


Harmony 4R144, 4R150

EN Instructions for use	5
DE Gebrauchsanweisung	14
FR Instructions d'utilisation	23
IT Istruzioni per l'uso	32
ES Instrucciones de uso	41
SV Bruksanvisning	50
PT Manual de utilização	59
NL Gebruiksaanwijzing	68
EL Οδηγίες χρήσης	77
JA 取扱説明書	86
ZH 使用说明书	95




Table of Contents

1.0 Application.....	6
1.1 Safety Considerations	6
1.2 Patient Information.....	6
2.0 Components	6
2.1 Replacement Parts and Accessories.....	7
3.0 Description	7
4.0 Installing the Harmony 4R134, 4R135, 4R144 and 4R150.....	8
4.1 Installing the Barbed Socket Connector.....	8
4.2 Attaching Socket and Adapter to Harmony.....	8
4.3 Attaching Pylon to Harmony	8
4.4 Connecting the Vacuum Line.....	8
5.0 Installing the Harmony 4R133	8
6.0 Adjusting the Harmony 4R144 and 4R150 Elastomer Rod Shock Absorber.....	9
7.0 Adjusting the Air Shock in the 4R133, 4R134, 4R1359.....	9
8.0 Monitoring the Elevated Vacuum Environment.....	10
8.1 Verifying the System if Vacuum Loss Occurs & Flow Chart	10
9.0 Technical Data	12
10.0 Warranty	12

1.0 Application

The Harmony Pump is a component of the Harmony Elevated Vacuum System, which also includes a Total Surface Weight Bearing Socket, Polyurethane Liner, Air Wick, and Sealing Sleeve. The Harmony is to be used solely for the prosthetic fitting of the lower limbs.

For Patient Weights of 50 to 100 kg. (110 to 220 lbs.)

- Harmony 4R133 Pump for use in some prosthetic feet (Fig. 5)
- Harmony 4R134 34mm Tube Clamp Adapter (Fig. 3)
- Harmony 4R135 Female Pyramid Receiver Adapter (Fig. 4)
- Harmony 4R144 Integrated Male Pyramid Connector and 30 mm Distal Tube Clamp Adapter (Fig. 1)

For Patient Weights of 100 to 150 kg. (220 – 330 lb.s)

- Harmony 4R150 Female Pyramid Receiver Adapter (Fig. 2)

1.1 Safety Considerations

- The Harmony Pump should be examined for proper function during routine check-ups by qualified personnel.
- The Harmony Pump may only be opened and repaired by an authorized Otto Bock Technician.
- The patient should be properly instructed in the use of the Harmony Pump and Harmony System as outlined in the Patient Information Section.

1.2 Patient Information

- Handle the Harmony Pump carefully and examine often for any signs of damage.
- Work with your prosthetist to make sure you understand how the Harmony System feels when functioning properly, when loss of vacuum occurs, and what the overall effects of applying an elevated vacuum to the limb.
- Read section 8.0 of this manual to understand how to perform basic system checks in case of loss of the vacuum.

<p>IMPORTANT: Proper maintenance and care require flushing 1 to 2 ounces of distilled water (about 5 to 10 pumps) through the Harmony® system after contact with salt water or perspiration.</p>

2.0 Components

The Harmony 4R144 and 4R150 pumps are delivered with the following components

1. Harmony Pump
2. Vacuum Gauge Assembly
3. Vacuum Line, Neoprene Boot, (2) Tube Collars and (1) Straight- and (1) 90°-Barbed Stainless Steel Socket Connectors
4. Instructions for Use SL=2900-4005

The Harmony 4R133, 4R134, 4R135 pumps are delivered with the following Components

1. Harmony Pump
2. Hand Pump
3. Tube Collar and (1) Straight- and (1) 90°-Barbed Stainless Steel Socket Connectors
4. Instructions for Use SL=2900-4005

2.1 Replacement Parts and Accessories

Harmony components can be ordered individually from Otto Bock using the following information.

Pc.	Description	Otto Bock TEC Part No.
1	Air Spring Pump Kit (for all units)	4Y352
1	Vacuum Line, 5 feet (for all units)	4Y309
10pk	Tube Collars (for all units)	4Y310
1	Tube Clamp Bar (for 4R134)	4Y323
1	Tube Clamp Screw, M5 x 22 Grade 12.9 (for 4R144)	501Z2=M5X22
1	34-30 mm reducing Split spacer (for 4R134)	4Y324
1	Valve Kit (for 4R134, 4R135)	4R138=5
5pk	Valve Kit, rubber air valves (for 4R144, 4R150)	4Y351-5
1	Elastomer Rod, 90 shore A, yellow (for 4R150)	4Y348
1	Elastomer Rod, 80 shore A red (for 4R144)	4Y347
1	Intake Valve with screen (for 4R144, 4R150) (Fig. 1, Pos.1)	4Y346
1	Expulsion Valve (for 4R144, 4R150) (Fig.1, Pos.2)	4Y345
1	1/8" Barbed Socket Connector/Filter (all units)	4Y344
1	1/8" Barbed Socket Connector 90° Elbow Connector/Filter (all units)	4Y350
3pk	Harmony Protective Boot (for 4R144, 4R150)	4Y319=3
3pk	Harmony Protective Boot (for 4R134, 4R135)	4Y319
3pk	Harmony Protective Boot (for 4R133)	4Y320
1	Harmony Complete Installation Kit (4R133, 4R134, 4R135)	4R130
1	Harmony Testing Kit, 4R133, 4R134, 4R135	4R131
1	Harmony 4R144, 4R150 Service Kit	

3.0 Description

The Otto Bock Harmony Pump combines the function of a vacuum pump and the function of a damping shock absorber in one unit. The 4R144 and 4R150 include the additional features of axial rotation and torsional shock absorption.

Function of the Vacuum Pump

The Harmony Pump is designed to generate an elevated vacuum to increase the linkage between user and prosthesis, and requires a sealed environment between the user's leg/liner, the socket, and the sealing sleeve. Once the environment is sealed, the Harmony Pump, which is compressed and released during normal walking, draws air from the socket through one-way valves and evacuates it through the pump. After multiple cycles of compression an elevated vacuum of 15-25 in. Hg. (508-847 millibars) is created in the sealed socket environment.

The Harmony Elevated Vacuum Socket System:

- Improves fluid exchange in the residual limb and helps to control daily volume fluctuations
- Increases proprioception
- Increases suspension to improve gait and overall function
- Reduces forces to the residual limb and liner
- Reduces moisture buildup on the residual limb

Function of the Shock Absorber

The Harmony shock absorption functions helps to reduce vertical and torsional forces in the extremities during walking. The result is an improvement for both the function and loading situation of the prosthesis wearer.

Vertical and Rotational Shock Absorption in the 4R144 and 4R150

Vertical shock absorption can be adjusted by the prosthetist by increasing or decreasing the compression of the elastomer rod. Rotational shock absorption is made possible by elastomer plates, which are not field serviceable. Vertical Shock Absorption in the 4R133, 4R134 and 4R135 can be adjusted by increasing or decreasing the amount of air pressure in the unit.

4.0 Installing the Harmony 4R134, 4R135, 4R144 and 4R150

Note: Failure to use the Harmony with a properly fitting Harmony socket system can compromise the limb, liner, and the Harmony pump.

4.1 Installing the Barbed Socket Connector

Note: To fabricate an airtight, laminated socket, blister-form a 3 mm PETG plastic sheet (575=1) as the first layer in the definitive socket lay-up.

To ensure an airtight system you must create an airtight channel between the Barbed Socket Connector and the interior of the socket. We recommend the following technique.

Prepare a channel in the socket by drilling a 3/8" (10 mm) hole in either the distal end or the posterior socket wall. Fill with five minute epoxy. Once epoxy has hardened, drill a 7/32 (5.5 mm) hole through the middle of the epoxy (Fig. 6). Use the Barbed Socket Connector to tap the hole, then remove and apply 5 minute epoxy to the threads, reinsert and hand tighten.

4.2 Attaching Socket and Adapter to Harmony

Pyramid Adapter style: Harmony 4R144

Tighten the set screws of the receiver onto pyramid and torque to 15NM (or to the manufacturer's specification if not an Otto Bock component). Use Loctite® 242 on screws for final adjustment.

4-Bolt Pattern Style: Harmony 4R134, 4R135 and 4R150

The Harmony 4R134, 4R135 and 4R150 pumps have a 4-bolt pattern on top for connecting the appropriate male or female adapters to the socket (i.e. 4R77, 4R23 or 4R51). Use the provided Cap Head Screws and connect the adapter to the Harmony. For final adjustment, apply Loctite® 242 to screws and torque to 115 in-lbs (13NM).

4.3 Attaching Pylon to Harmony

30mm Tube Clamp (4R144): Insert the 2R2 pylon into tube completely. Torque the tube clamp screw to 88 in-lbs (10 Nm) for the final adjustment.

Female Pyramid Receiver Style (4R135 and 4R150): Tighten the set screws of the pyramid receiver onto the attachment and torque to 133 in-lbs (15Nm.). Use Loctite® 242 on screws for final adjustment.

34mm Tube Clamp (4R134): Insert 2R76 pylon into tube clamp completely. Torque the tube clamp screw to 44 in-lbs (5 Nm) for the final adjustment.

4.4 Connecting the Vacuum Line

To connect the Vacuum Line to the socket, cut line to length proximally, moisten end of line, slide Retention Collar over Line, attach Line to Barbed Socket Connector, and slide Retention Collar into place (Fig. 7).

5.0 Installing the Harmony 4R133

- 5.1** To remove original shock, remove the cap head screw from the housing of the foot with a 4mm Allen wrench.
- 5.2** Push the pin through the shock and remove shock.
- 5.3** Repeat steps 5.1 and 5.2 for the other end of the shock.
- 5.4** Remove the original shock from the housing of the foot.
- 5.5** Remove the washers from the original shock.
- 5.6** Remove the bushings from the original shock and insert bushings in the Harmony 4R133.

- 5.7 Re-use the original washers on the proximal end of the Harmony 4R133.
- 5.8 Place provided washers supplied with Harmony 4R133 on both sides of the distal bushings.
- 5.9 Install the Harmony 4R133 in the foot housing.
- 5.10 Push the pins through the bushings.
- 5.11 Reinstall the cap head screws in the pin and tighten with a 4 mm Allen wrench.
- 5.12 To connect the Vacuum Line to the socket, see section 4.4.

6.0 Adjusting the Harmony 4R144 and 4R150 Elastomer Rod Shock Absorber – PRACTITIONERS ONLY

The preload on the elastomer rod in these units can be adjusted to accommodate different patient weights. This function should be adjusted for the patient's normal gait. In order for the vacuum pump to function properly, the Harmony 4R144 and 4R150 must be adjusted so that the shaft extends completely when unweighted (with no air pressure in the air chamber), and so that it compresses between 5 and 7 mm when weighted.

To replace the Elastomer Rod:

1. Locate Blue Cup inside shaft (Fig. 8). Screw in completely using a 3/8" (10mm) Allen wrench
2. To set Elastomer Rod tension, release pressure on the Rod by backing out the Blue Cup. See suggested starting points below. Because the resistance of the elastomer rods can change slightly over time, we recommend it be readjusted during routine checkups.

4R144 only (Red Elastomer Rod)

Patient Weight (lbs / kgs)	# of turns out (counterclockwise) from fully inset position
120 lb. / 50kg.	4.5
140 lb. / 60 kg.	4.0
160 lb. / 70 kg.	3.5
180 lb. / 80 kg.	3.0
200 lb. / 90 kg.	2.5
220 lb. / 100 kg	2.0

4R150 only (Yellow Elastomer Rod)

Patient Weight (lbs / kgs)	# of turns out (counterclockwise) from fully inset position
220 lb. / 100kg.	4.5
240 lb. / 110 kg.	4.0
260 lb. / 120 kg.	3.5
280 lb. / 130 kg.	3.0
300 lb. / 140 kg.	2.5
320 lb. / 150 kg.	2.0

Important: Do not unscrew more than 5 turns, as this may prevent the shaft from extending completely and compromise the vacuum pump function.

7.0 Adjusting the Air Shock in the 4R133, 4R134, 4R135 – USERS AND PRACTITIONERS

Resistance in the Air Shock can be modified with a hand pump when high-level activities are encountered (Fig. 9). Please be aware that increased resistance in the Air Shock may limit the travel of the pump piston, requiring more compression cycles for you to attain the desired elevated vacuum during normal walking.

Note: When adjusting pressure in the Air Shock, we recommend that you begin at 0 psi as a starting point. Please see section 9.0 for recommended air pressure.

8.0 Monitoring the Elevated Vacuum Environment

A Vacuum Gauge Assembly is included with the 4R144 and 4R150 Pumps.

We recommend purchasing the 4R130 Harmony Complete Installation Kit, or the 4R131 Harmony Testing Kit for monitoring and repair needs for the 4R133, 4R134, 4R135 units.

To verify that the Harmony System can maintain an elevated vacuum environment, it is necessary to temporarily install a vacuum gauge. To do this, remove the Vacuum Line from the Harmony Intake Valve (Fig. 10). Next, take the Vacuum Gauge Assembly and attach the Vacuum Line (short side) to the Inlet Valve (Fig.1, Pos. 1) and connect the socket-side Vacuum Line to the open barb on Vacuum Gauge Assembly (Fig.11).

While the patient is standing and wearing the properly fitting Harmony System, cycle the Harmony Pump to generate a vacuum between 15-25 in. Hg. (508-847 millibars). Stop cycling the pump and verify that the vacuum is maintained when it is not being cycled. If the vacuum is maintained, the system is functioning correctly. If vacuum loss does occur, it will be necessary to check the system as shown in Section 8.1.


8.1 Verifying the System if Vacuum Loss Occurs – See Flow-Chart

To determine if the vacuum loss exists in Zone 1, 2 or 3 (Fig. 12) install gauge as described in section 8.0 (Fig. 11).

Cover the hole on the inside of the socket with a piece of tape that is not air-permeable (such as PVC tape). Cycle the pump to reach 15-25 in. Hg. (508-847 millibars). Check to see if vacuum level is maintained when the pump is not being cycled.

Continue to the Flow-Chart

FLOW CHART

Note: To clean the Barbed Socket Connector: Disconnect the tubing from the Barbed Socket Connector and gently blow air through the fitting back into the socket.

9.0 Technical Data

Harmony-4R133

Weight: .68 lbs/309 grams

Weight Limit: 220 lbs/100 kgs

Air spring operating psi: 40-120 psi/1.35-4.06 bar

Material: Aluminum

Harmony-4R134

Weight: 1.6 lbs/726 grams

Height: 5 1/8 in/130 mm

Diameter: 3.2 in/81mm

Weight Range: 110-220 lbs/50-100 kgs

Air spring operating psi: 0-20 psi/0-.68 bar

Tube clamp screw: 44 in/lbs (5 Nm)

Adapter 4-Hole pattern: 120 in/lbs (13 Nm)

Material: Aluminum/Titanium

Harmony-4R135

Weight: 1.6 lbs/726 grams

Height: 4 5/8 in/117 mm

Diameter: 3.2 in/81mm

Weight Range: 110-220 lbs/50-100 kgs

Air spring operating psi: 0-20 psi/0-.68 bar

Receiver Set Screw Torque: 133 in/lbs (15 Nm)

Adapter 4-Hole pattern: 120 in/lbs (13 Nm)

Material: Aluminum/Titanium

Harmony-4R144

Weight: 1.1 lbs/520 grams

Height: 5 1/8 in/130.6 mm

Weight Range: 110-220 lbs/50-100 kgs

Tube Clamp Screw Torque: 88 in/lbs (10 Nm)

Air spring operating psi: 0-20 psi/0-.68 bar

Material: Aluminum/Stainless Steel

Harmony-4R150

Weight: 1.4 lbs/640 grams

Height: 5.0 in/127 mm


Weight Range: 220-330 lbs/100-150 kgs

Receiver Set Screw: 133 in/lbs (15 Nm)

Adapter 4-Hole pattern: 120 in/lbs (13 Nm)

Air spring operating pressure: 0-40 psi/0-1.36 bar

Material: Aluminum/Stainless Steel

10.00 LIMITED WARRANTY

Otto Bock HealthCare (Otto Bock) warrants all of its products, to the original purchaser, to be free from defects in materials and workmanship. This warranty applies, subject to normal wear and tear, when the products are used as intended, without unapproved modifications, following all Otto Bock instructions and requirements; and when they are fitted by or under the direct supervision of certified/licensed practitioners who meet all Otto Bock product-specific training requirements for a given product.

The duration of our Limited Warranties varies by product, but all are effective from the date of fitting to a patient.

The duration of the Limited Warranty for some components may be doubled when incorporated into a finished device in which all of the components are manufactured by Otto Bock and installed in conformance with all Otto Bock recommendations. Otto Bock's sole obligation under this Limited Warranty shall be to repair or replace at no charge or to refund the cost of the system or component thereof to the original purchaser, at Otto Bock's sole discretion.

The Limited Warranty for the Harmony P2 is 12 months, and qualifies for extension to 24 months when in finished devices utilizing all Otto Bock components.

THE EXPRESS WARRANTIES SET FORTH HEREIN ABOVE ARE IN LIEU OF ALL OTHER WARRANTIES EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, AND ALL SUCH OTHER WARRANTIES ARE HEREBY DISCLAIMED AND EXCLUDED BY OTTO BOCK.

IN NO EVENT SHALL OTTO BOCK'S LIABILITY OF ANY KIND INCLUDE ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, EVEN IF OTTO BOCK SHALL HAVE BEEN ADVISED OF THE POSSIBILITY OF SUCH POTENTIAL LOSS OR DAMAGE.

Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation may not apply to you.

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Otto Bock HealthCare GmbH