
## 7E7




<b>DE</b> Gebrauchsanweisung .....	5
<b>EN</b> Instructions for use .....	12
<b>FR</b> Instructions d'utilisation .....	19
<b>IT</b> Istruzioni per l'uso .....	26
<b>ES</b> Instrucciones de uso .....	33
<b>PT</b> Manual de utilização .....	41
<b>NL</b> Gebruiksaanwijzing .....	48
<b>SV</b> Bruksanvisning .....	55
<b>PL</b> Instrukcja użytkowania .....	62
<b>HU</b> Használati utasítás .....	69
<b>CS</b> Návod k použití .....	77
<b>HR</b> Upute za uporabu .....	84
<b>EL</b> Οδηγίες χρήσης .....	91
<b>RU</b> Руководство по применению .....	98






● Einzelteile-Pack 7D2  
 7D2 Single-Component Pack  
 Kit de pièces de rechange 7D2  
 Imballo dei singoli componenti 7D2  
 7D2 Kit componentes  
 Service-Set 7D2  
 7D2 Onderdelenpakket  
 Conjunto de peças de reposo 7D2  
 Alkatrész-csomag 7D2-hez  
 Sada jednotlivých dílů 7D2  
 Комплект деталей 7D2

▲ Mindestmenge  
 Minimum order quantity  
 Quantité minimum  
 Quantità minima  
 Cantidad mínima  
 Minimikvantitet  
 Minimum aantal  
 Quantidade minima  
 Minimum mennyiség  
 Minimální odběr, množství  
 Минимальное количество

■ Einzelteile  
 Single components  
 Pièces à l'unité  
 Singoli componenti  
 Componentes  
 Separata delar  
 Onderdelen  
 Componentes  
 Alkatrészek  
 Jednotlivé díly  
 Отдельные детали

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<b>⚠ WARNING</b>	Warnings regarding possible risks of severe accident or injury.
<b>⚠ CAUTION</b>	Warnings regarding possible risks of accident or injury.
<b>NOTICE</b>	Warnings regarding possible technical damage.
<b>INFORMATION</b>	Additional information on the fitting / use.

## INFORMATION

Last update: 2021-02-17

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document in a safe place.

## 1 Single components (Fig. 1)

### 1.1 Components

- |   |   |
|---|---|
| (1) 2R25=16-400 Aluminium Tube              | (17) 4Z30 Safety Plate                              |
| (2) 4G80=N4 Upper Joint Section – Titanium– | (18) 501S27=M4x5 Oval Countersunk Head Screw        |
| (3) 4G81 Lower Joint Section – Titanium–    | (19) 7Z55 Stop                                      |
| (4) 7Z42=C Stop Wedge                       | (20) 4B130 Collar Bushing (2 pieces)                |
| (5) 7Z43 Adjustment Wedge                   | (21) 4G194=N Bottom Part                            |
| (6) 7Z49 Guide Pin                          | (22) 501F9=M8x30 Six Lobe Pan Head Screw (2 pieces) |
| (7) 7Z50=H Tappet                           | (23) 501Z2=M5x30 Cap Screw DIN 912                  |
| (8) 7Z51 Guide Bushing                      | (24) 4A54=4x15.8 Axis                               |
| (9) 7Z52 Adjustment Nut                     | (25) 7Z60 Lock Pin                                  |
| (10) 513D36 Extension Assist Spring         | (26) 7Z58=N Lamination Dummy                        |
| (11) 7Z53 Lamination Plate                  | (27) 501T39=M4x12 Cap Screw                         |
| (12) 7Z54 Clamping Jaw (2 pieces)           | (28) 4Z48=1 Brass Washer (not illustrated)          |
| (13) 506G1=M8x10 Set Screw                  | (29) 709Z11 Bit T40 (Torx 40) (not illustrated)     |
| (14) 7Z57 Clamping Disc (2 pieces)          | (30) 4G507 Intermediate Piece                       |
| (15) 4A42=46 Axis                           |   |
| (16) 4A43 Axis Nut                          |   |

## 1.2 7D2 Single-Component pack for 7E7 (•)

Consisting of:


1 Tappet 7Z50=H; 1 Guide Bushing 7Z51; 1 Extension Assist Spring 513D36; 1 Safety Plate 4Z30; 1 Oval Countersunk Head Screw 501S27=M4x5; 1 Stop 7Z55; 2 Pan Head Screws 501F9=M8x30; 1 Cap Screw 501T39=M4x12; 1 Cap Screw 501T39=M4x14; 1 Cap Screw 501Z2=M5x30, 1 Lock Pin 7Z60

## 2 Description

### 2.1 Indications for use

The 7E7 Modular Hip Joint with Internal Extension Assist is to be used **exclusively** for the exo-prosthetic fitting of patients with hip amputations such as intertrochanteric amputation, hip disarticulation or hemipelvectomy.

### 2.2 Field of application

Field of application according to the **Ottobock MOBIS Mobility System**:

Recommended for amputees with Mobility Grades **2 and 3**

(restricted outdoor walkers and unrestricted outdoor walkers).

Approved for a **patient weight of up to 100 kg**.

### 2.3 Combination possibilities

#### CAUTION

##### **Failure to observe the manufacturer requirements regarding combination possibilities**

Injuries, malfunctions or product damage due to unallowable combination of prosthesis components

- Based on the instructions for use of all prosthetic components used, verify that they may be combined with each other and are approved for the patient's field of application.

#### INFORMATION

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.

#### **Recommended combinations**

Prosthetic knee joints | Kenevo 3C60, C-Leg 3C98-3, Genium 3B1-3, Genium X3 3B5-3

#### **Unallowable combinations**

Prosthetic knee joints | Kenevo 3C60=ST, C-Leg 3C88-3, Genium 3B1-3=ST, Genium X3 3B5-3=ST

### 2.4 Reuse and lifetime

#### CAUTION

##### **Reuse on another patient**

Fall due to loss of functionality as well as damage to the product

- Only use the product for a single patient.

## CAUTION

### Exceeding the lifetime

Fall due to change in or loss of functionality and damage to the product

- Ensure that the maximum lifetime defined in this section is not exceeded.

This product was tested by the manufacturer with 2 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 3 years.

## 2.5 General safety instructions

Please make your patients aware of the following safety instructions:

## WARNING

**Risk of accident when driving a motor vehicle.** The ability of prosthesis users to drive a vehicle is determined on a case-by-case basis. Criteria include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's individual abilities.

All persons are required to observe their country's national and state driving laws when operating motor vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorised test centre.

For maximum safety and convenience, Ottobock recommend that, at the very least, a specialist evaluate the need for any adaptations to the car (such as by installing a steering fork, automatic shift). Risk-free driving should be ensured even when the prosthesis is not functioning.

## CAUTION

**Risk of injury as a result of excessive strain.** The Ottobock prosthetic components described in these instructions for use have been developed for activities of daily living and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's safety.

Should prosthetic components be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the components inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Service.

## CAUTION

**Risk of pinching.** During hip joint flexion there is the risk of pinching between the socket and hip joint. It must be explicitly pointed out to the prosthesis wearer that there must not be any body parts or objects in this area (e.g. mobile phone in the trouser pocket). It is recommended to use a cosmetic cover for the joint with every fitting, even in case of trial fittings.

## NOTICE

**Damage due to corrosion.** Prosthetic components must not be exposed to corrosive elements such as fresh water, salt water and acids. Using this medical device in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

## NOTICE

**Damage due to incorrect environmental conditions.** Prosthetic components must not be subjected to intense smoke, dust, vibrations, shocks or high temperatures. Do not allow debris or liquids to get into the prosthesis and its components. Failure to follow these instructions can lead to malfunctioning and damage to the prosthesis.

## NOTICE

**Damage caused by improper cleaning agents.** Do not use any aggressive cleaning agents since they could cause damage to bearings, seals and plastic parts.

### 2.6 Function

The joint features small structural height (= laminate thickness), which helps to minimise 'pelvic tilt' while the patient is sitting. It has infinitely adjustable abduction/adduction, flexion/extension and rotation. The adjustment of the hip joint in relation to the knee axis is performed by rotating the aluminium tube.

The tension of the integrated extension assist is also infinitely adjustable. The extension assist spring can be accessed by removing the aluminium tube.

The lamination plate can be removed from the joint section and should only be adjusted while the tuber is in horizontal position.

## 3. Handling

### 3.1 Pelvic socket fabrication

#### CAUTION

**Risk of falling due to incorrect lamination of the lamination plate.** Proper lamination technique is necessary to guarantee stability and comfort for the user. It is particularly important that the lamination plate is sufficiently reinforced with carbon-fibre cloth and fibreglass matting.

#### CAUTION

**Risk of health impairment when processing fibreglass.** Fibreglass can irritate the eyes, respiratory tract and skin. Where appropriate, use a particle filter mask or exhaust equipment.

#### 3.1.1 Preparatory work

- Drill three evacuation ducts into the plaster positive to ensure intimate fit.
- Align the lamination plate on the plaster positive, shorten as necessary and round the edges.
- If there is large tissue loss on the residual limb, or in case of hemipelvectomy, first restore symmetry by using Pedilen® rigid foam.
- Pull over the form-weld inner bag and evacuate.
- Manufacture the external bag with pouring hole located on the side of the residual limb.
- Cut off the filling hoses for the hard/soft lamination technique.

#### 3.1.2 Reinforcement of the pelvic socket and lamination plate

The following reinforcement material shall be used depending on the type of loading expected:

- 623T3 Perlon Stockinette (10-12 layers)

- 616G12 Carbon Fibre Cloth (2 - 4 layers)
- 616G4 Fibreglass Mat (2 layers)

The following working steps are necessary:

- Cut a length of Perlon stockinette twice as long as the plaster positive, pull over the plaster positive, twist on the non-amputated side and reflect.
- Pull over 2 other Perlon layers.
- Reinforce the socket in the seat area with carbon fibre cloth and fibreglass mat.
- Close the threaded holes in the lamination plate with 636K8 Plastaband.
- Interweave 699B1 Fibreglass Roving through the lamination plate holes.
- Align the lamination plate to the marking on the plaster positive and apply fibreglass mat and carbon fibre cloth.
- Complete the reinforcement with Perlon stockinette while inserting the filling hoses, Trolene separation flap and threaded plates.
- Pull the external bag over the reinforced plaster positive and thread the filling hoses through the pouring hole.

### 3.1.3 Lamination technique

By pouring 617H17 Orthocryl Flexible Resin into the ventrally and dorsally inserted filling hoses, these sections can be formed as desired. Pour 617H19 Orthocryl Lamination Resin into the affected side.

It is important that the lamination and evacuation steps are performed separately and, if possible, by 2 prosthetists, to ensure that the resin remains in the desired position during the process

### 3.1.4 Completion and setting-up

- Once the resin has hardened, cut open the pelvic socket at its brim.
- Use a knife to cut open the outside and inside flaps along the Trolene brim.
- Remove the pelvic socket from the plaster positive and grind the brim smooth.
- Expose the threaded hole of the lamination plate and mount the hip joint to the pelvic socket.

### 3.1.5 Using the 7Z58=N lamination dummy

After having completed the trial fitting of the pelvic socket with the mounted hip joint, the hip joint must be disassembled again and a protective or connection cap be fabricated using the enclosed lamination dummy in order to ensure the function of the hip joint and for the later connection to a cosmetic foam cover. To position the lamination dummy/ protective cap properly, the protruding guide pins of the lamination dummy must be inserted into the two threaded holes provided in the lamination plate.

## 3.2 Alignment

Figure 2 shows the alignment recommendation for a pelvis prosthesis with 7E7 Hip Joint, 4R56 Angled Tube Clamp Adapter, 3R60 Knee Joint and 1C30 Trias.

The Angled Tube Clamp Adapters are designed especially for adjustable connection between the tube of the hip joint and pyramid adapter of the knee joint. The Angled Tube Clamp Adapter has a 10° angle (4R56), a 20° angle (4R56=1) or a 30° angle (4R56=2) depending on the position of the tube of the anteriorly located hip joint.

## 3.3 Adjustment options

The 7E7 Hip Joint is adjustable in three planes: abduction-adduction, flexion-extension and hip axis rotation (Fig. 3).

### 3.3.1 Abduction - adduction to the pelvic socket

Using a 710D4 Torque Wrench and 709Z11 Bit T40 (29), loosen the two 501F9=M8x25 Pan Head Screws (22) (Fig. 4) until the hip joint can move with respect to the pelvic socket. Set the required position  $\pm 10^\circ$  (Fig. 5.1).

### 3.3.2 Rotation position of the hip axis

The rotation adjustment can be set between the Upper Joint Section (2) and Bottom Plate (21). To do this, loosen the two pan head screws with the 710D4 Torque Wrench and 709Z11 Bit T40 (29) (Fig. 4) and additionally loosen the Cap Screw (27) with a 3-mm Allen wrench (Fig. 6). Setting range  $\pm 5^\circ$  (Fig. 5.2, 5.3).

#### INFORMATION

Use the Cap Screw (27) only for adjustment of the  $\pm 5^\circ$  rotation with respect to the pelvic socket. Maximum torque: **M=2 Nm** (Fig. 5.2 + 5.3 and 6).

For trial walking or upon completion of the prosthesis, tighten both 501F9=M8x30 Pan Head Screws with a torque of **M=25 Nm** (Fig. 4). Secure the Pan Head Screws with 636K13 Loctite while completing the prosthesis.

### 3.3.3 Flexion – extension

The desired flexion-extension setting can be attained with the Cap Screw (23) with a 4-mm Allen wrench (Fig. 7). When you are finished with the adjustment, insert the Lock Pin (25) (Fig. 9, arrow).

### 3.3.4 Adjusting the tension of the hip extension assist

Loosen the Set Screw on the Lower Joint Section (Fig. 8). Then pull out the Aluminium Tube (1). Rotate the Adjustment Nut (9) to increase or decrease the spring tension (Fig. 9).

Then insert the Aluminium Tube (1) into the Lower Joint Section all the way to the stop and clamp in position by tightening the Set Screw (13). **Use 710D4 Torque Wrench. Torque: 15 Nm. Secure the set screw with 630K13 Loctite.**

#### CAUTION

**Risk of falling when using the tube clamp for length adjustment.** The minimum clamp length of the tube must be observed. Make sure that the tube is fully inserted into the tube clamp.

#### INFORMATION

Only the 2R25=16-400 Aluminium Tube may be used. Please note the different article number when purchasing a replacement tube later on: 2R30.

Use Lamination Dummy (26) to fabricate a connection cap for the cosmetic foam cover.

## 3.4 Service

We recommended that you check the Modular Hip Joint after several weeks and adjust if necessary.

### 3.4.1 Replacing the bearing bushings:

- Press out worn bearing bushings.
- Press in two 4B130 Bushings at right angles with the side face of 4G80=N4.
- Ream the bushings using 727S9=12G7 Reamer.
- Lubricate the bushings with 633F30 Special Lubricant.

### 3.4.2 Lubrication

For lubricating the plastic components 633G6 Lubricant is recommended.

For lubricating metal parts sliding against each other such as Tappet (7) and Bearing Axle (24) use 633F16 Gliding Paste. For lubricating 4B130 Collar Bushing (20) and Axis (15) use 633F30.

### 3.5 Technical information

Distal connection	Tube clamp
Max. flexion angle	140°
Weight of the joint	340 g
Weight of connection technology	280 g
System height	33 mm
Material	Titanium
Mobility Grade	2 + 3
Max. body weight	100 kg

## 4 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

## 5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### 5.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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