



ottobock.

CE

9X50/9X51

DE Gebrauchsanweisung (Fachpersonal)	3
EN Instructions for use (qualified personnel)	11
FR Instructions d'utilisation (Personnel spécialisé)	20
IT Istruzioni per l'uso (Personale tecnico specializzato)	29
ES Instrucciones de uso (Personal técnico especializado)	38
PT Manual de utilização (Pessoal técnico)	47
NL Gebruiksaanwijzing (Vakmensen)	56
SV Bruksanvisning (Fackpersonal)	65
DA Brugsanvisning (Faguddannet personale)	73
EL Οδηγίες χρήσης (Τεχνικό προσωπικό)	81

INFORMATION

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INFORMATION

Date of last update: 2020-12-23

- ▶ Please read this document carefully before using the product and observe the safety notices.

- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product description

2.1 Function

9X50 linear control element function

The 9X50 linear control element is installed in the definitive socket. The switching points are reached by pulling on the bow. When the 9X50 linear control element is not under load, the control signal is 0 volts. When the bow is pulled to maximum, the control signal corresponds to the maximum electrode signal of an electrode.

9X51 4-stage control element function

Pulling on the bow of the 9X51 4-stage control element generates a control signal. Output signals of various strengths are simulated depending on the ratchet position, making it possible to switch between the various components with corresponding configuration of the prosthesis.

2.2 Combination possibilities

This product can be combined with the following Ottobock components:

- DynamicArm: 12K100N=* / 12K110N=*
- ErgoArm Hybrid plus: 12K44=*
- ErgoArm Electronic plus: 12K50=*
- AxonMaster 13E500
- System Electric Hand DMC plus: 8E38=6, 8E39=6
- System Electric Hand Digital Twin: 8E38=7, 8E39=7
- Sensor Hand Speed: 8E38=8, 8E39=8
- MyoHand VariPlus Speed: 8E38=9, 8E39=9, 8E41=9
- Transcarpal Hand DMC plus: 8E44=6
- System Electric Greifer DMC VariPlus: 8E33=9, 8E34=9
- 8E70=* bebionic hand EQD

- 8E71=* bebionic hand Short Wrist
- 8E72=* bebionic hand Flex
- 10S17 electric wrist rotator
- 13E205 MyoRotronic

3 Intended use

3.1 Indications for use

The product is to be used solely for upper limb prosthetic fittings.

3.2 Conditions of use

The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

Permissible ambient conditions are described in the technical data (see page 19).

3.3 Indications

- For unilateral or bilateral amputation
- Dysmelia of the forearm or upper arm

3.4 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

3.5 Qualification

The fitting of a patient with the product may only be carried out by O&P professionals who have been authorised with the corresponding Ottobock training.

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
NOTICE	Warning regarding possible technical damage.

4.2 Structure of the safety instructions

WARNING

The heading describes the source and/or the type of hazard

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in the event of failure to observe the hazard
- > E.g.: Consequence 2 in the event of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

CAUTION

Non-observance of safety instructions

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety instructions and the stated precautions in this accompanying document.

CAUTION

Use of a damaged product

Injury due to loss of product functionality.

- ▶ Before use, perform a visual inspection to check if all parts of the product are undamaged.
- ▶ Replace the product immediately in case of damage.

CAUTION

Independent manipulation of the product

Injury due to product damage or malfunction.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ Opening and repairing the product or repairing damaged components is not permitted.

CAUTION

Signs of wear on the product

Injury due to faulty control or malfunction of the product.

- ▶ If a malfunction of the product is noted, it has to be replaced.
- ▶ Opening and repairing the product or repairing damaged components is not permitted. .

⚠ CAUTION

Faulty adjustment

Injury due to faulty control or malfunction of the product.

- ▶ When mounted into a provisional socket, the product may only be worn within an orthopaedic workshop or as an inpatient in a rehabilitation facility.

⚠ CAUTION

Penetration of liquids

Injury due to faulty control or malfunction of the product.

- ▶ Make sure that no liquids can penetrate into the product.
- ▶ Do not use the product for bathing prostheses.

⚠ CAUTION

Damage to the connection cable

Injury due to faulty control or malfunction of the product.

- ▶ When installing the connection cable, make sure there are no tight angles or kinks.
- ▶ Replace damaged connection cables immediately.

⚠ CAUTION

Faulty adjustment

Injury due to faulty control or malfunction of the product.

- ▶ A load may only be applied to the product in the pulling direction.
- ▶ Please note that the cable travel of the product is limited. To prevent excessive strain on the product, make sure the harness is correctly fitted and accommodates the patient's movements.
- ▶ In the course of service, check the bow and housing of the product for damage.

4.4 Patient information

WARNING

Non-observance of safety notices

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

WARNING

Operating vehicles and machinery

Injury due to faulty control or malfunction of the product.

- ▶ Driving vehicles of any kind or operating machinery is strictly prohibited while wearing a provisional socket and the product.

CAUTION

Penetration of dirt and humidity into the product

Injury due to faulty control or malfunction of the product.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

CAUTION

Independent user changes or modifications made to the product

- > Injury due to faulty control or malfunction of the product.
- > Destruction of the product.
- ▶ Only have any changes or modifications to the product carried out by a certified Ottobock Service Centre.

CAUTION

Mechanical loads

Injury due to faulty control or malfunction of the product.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product and its housing for visible damage before each use (e.g., cracking or breakage). Should the product be damaged or not function properly, please contact your orthopaedic technician immediately.
- ▶ When the product encounters a noticeable stop, do not continue increasing the pull force.

⚠ CAUTION**Remaining in areas outside the allowable temperature range**

Injury due to faulty control or malfunction of the product.

- ▶ Avoid remaining in areas outside the allowable temperature range (see page 19).

⚠ CAUTION**Signs of wear on the product**

Injury due to faulty control or malfunction of the product.

- ▶ Since signs of wear can occur with all moving parts, annual service is required.

NOTICE**Improper care**

Damage or breakage through the use of acetone, petrol or similar solvents.

- ▶ Clean the product and the contacts with a damp cloth and mild soap only.

5 Scope of Delivery and Accessories

5.1 Scope of delivery

- 1 pc. 9X50 linear control element
- 1 pc. 9X51 4-stage control element
- 1 pc. Instructions for use (qualified personnel)

5.2 Accessories

The following components are not included in the scope of delivery and must be ordered separately:

- 13E129=* electrode cable

6 Preparing the product for use

6.1 Required tools

- 633F11 silicone grease

6.2 Connecting the electrode cable

- 1) Shorten the electrode cable as needed.
- 2) Push the electrode cable with the grey side up into the plug connection all the way to the stop (see fig. 1).

- 3) Apply sufficient silicone grease to the inside of the control element housing and the contact pins (see fig. 2).
- 4) Connect the plug connection between the housing and electrode cable plug.
- 5) Push the plug on the electrode cable as deep as possible into the product using a screwdriver (see fig. 3).
- 6) Remove any silicone grease that seeps out.

Electrode cable replacement

- 1) Pull off the electrode cable at right angles to the product (see fig. 4).
- 2) Cut off at least 5 mm of the electrode cable and push it back into the plug connection.

7 Use

The 9X50 linear control element and the 9X51 4-stage control element are attached between the shoulder harness and definitive socket (see fig. 5).

8 Legal information

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

8.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

8.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

9 Technical data

Ambient conditions	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing

General information	
Reference number	9X50/9X51
Product service life	5 years
Cable travel	8 mm
Max. actuating force	10 N
Output voltage	0 - 2.5 V
Power supply	757B20, 757B21, 757B35=*, 757B15

10 Appendices

10.1 Symbols Used

Declaration of conformity according to the applicable European directives

Lot number (PPPP YYYY WW)

PPPP – plant

YYYY – year of manufacture

WW – week of manufacture

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

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