3R106-Pro, 3R106-Pro=KD, 3R106-Pro=ST

DE Gebrauchsanweisung	7
EN Instructions for use	20
FR Instructions d'utilisation	32
IT Istruzioni per l'uso	46
ES Instrucciones de uso	59
PT Manual de utilização	73
NL Gebruiksaanwijzing	87
SV Bruksanvisning	100
DA Brugsanvisning	112
NO Bruksanvisning	125
FI Käyttöohje	138
PL Instrukcja użytkowania	150
HU Használati utasítás	164
CS Návod k použití	176
HR Upute za uporabu	189
TR Kullanma talimatı	201
RU Руководство по применению	214
JA 取扱説明書	228
ZH 使用说明书	240

1

①

②

③



④

2Z11=KIT

⑤

⑥

⑦

⑧

⑨





⑩




⑪

INFORMATION

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INFORMATION

Date of last update: 2021-01-11

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

1.1 Construction and Function

The 3R106-Pro* Prosthetic Knee Joint is a polycentric prosthetic knee joint. It is available in the following versions:

- 3R106-Pro (with proximal pyramid adapter; in combination with the 4R156* Tube Clamp Adapter for hip disarticulation)
- 3R106-Pro=KD (with proximal lamination anchor for knee disarticulation)
- 3R106-Pro=ST (with proximal threaded adapter for long residual limb prostheses)

The prosthetic knee joint features a large flexion angle (175°) (see fig. 3). It is equipped with auto-adaptive, pneumatic swing phase control for walking at different speeds. Flexion ① and extension damping ② can be adjusted separately and individually (see fig. 4).

The upper joint section ① and lower joint section ② are connected by the two anterior links ③ and the posterior link ④ to form a four-bar linkage. In extension, the instantaneous centre of rotation ⑤ is located above the joint and behind the weight bearing line, thus providing knee stability in the stance phase (see fig. 2).

The distal portion of the pneumatic cylinder is fabricated as a tube clamp and accommodates the provided tube adapter.

The integrated extension assist spring is part of the swing phase control mechanism and ensures extension of the prosthetic knee joint. If required, it can be replaced by the included weaker spring. Here it is important to ensure that the patient always achieves full extension while walking.

1.2 Combination possibilities

INFORMATION

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.

Prosthetic components that are especially well suited for combination with the product are listed below.

Designation	Reference number
Tube Clamp Adapter	4R156* (for hip disarticulation)
Foam Cover	3S107 (for 3R106-Pro and 3R106-Pro=ST, if needed also 3R106-Pro=KD) 3R27* (for hip disarticulation) 6R6 (for 3R106-Pro=KD)

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

3R106-Pro*

The product is recommended for mobility grade 2 (restricted outdoor walker) and mobility grade 3 (unrestricted outdoor walker). Approved for a body weight of **up to 125 kg (275 lbs)**.

2R49

Approved for a body weight of up to **125 kg**.

2.3 Environmental conditions

Allowable environmental conditions

Temperature range for use: -10 °C to +60°C

Allowable relative humidity 0 % to 90 %, non-condensing

Unallowable environmental conditions

Mechanical vibrations or impacts

Perspiration, urine, fresh water, salt water, acids

Dust, sand, highly hygroscopic particles (e.g. talcum)

2.4 Lifetime

This product was tested by the manufacturer with 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 5 years.

3 Safety

3.1 Explanation of warning symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

3.2 General safety instructions

CAUTION

Reuse on another patient

Risk of injury due to loss of functionality as well as damage to the product

► Only use the product for a single patient.

CAUTION

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

► Use the prosthetic components according to the MOBIS classification.

CAUTION

Improper combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

⚠ CAUTION

Use under restricted environmental conditions

Risk of injury due to damage to the product

- ▶ Do not expose the product to restricted environmental conditions.
- ▶ If the product has been exposed to restricted environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

⚠ CAUTION

Exceeding the service life

Risk of injury due to change in or loss of functionality and damage to the product

- ▶ Ensure that the approved service life is not exceeded.

⚠ CAUTION

Reaching into the area of the joint mechanism

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.

⚠ CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Noticeable changes in functionality may include poor response, incomplete extension, diminished swing phase control or stance phase stability, abnormal noises, etc.

4 Scope of delivery

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

Fig.		Item no.	Quantity	Designation	Reference number
1	-	①	1	Prosthetic Knee Joint	-
1	■	②	1	Tube Adapter	2R49
1	■	③	0	Adjustment Wrench	710H10=2X3
-	■		1	Compression Spring, heavy duty	513D83=1.4X12.6X66
1	■	④	1	Compression Spring, light duty	513D83=1.1X12.9X66

Fig.		Item no.	Quantity	Designation	Reference number
–	■		1	Instructions for Use	647G208=all_INT
The following parts are only for:				Prosthetic Knee Joint	3R106-Pro
1	■	⑤	1	Pyramid Dome Protektor Kit	2Z11=KIT
The following parts are only for:				Prosthetic Knee Joint	3R106-Pro=KD
1	■	⑥	1	Lamination Anchor	4G70
The following parts are for reorder only (not included in the scope of delivery!)					
1	●		0	Single Component Pack	4D3
	Consisting of:				
	■	⑦	1	Bumper	4Z81=4X12
	■	⑧	4	Truss Head Screw	501T1=M5X16
	■	⑨	4	Two-Hole Nuts	502R1=M5X16
	■	⑩	2	Set Screw	506G3=M8X12-V
	■	⑪	1	Set Screw	506G3=M8X10

5 Preparing the product for use

5.1 Information on fabrication of a prosthesis

CAUTION

Incorrect alignment, assembly or adjustment

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- Observe the alignment, assembly and adjustment instructions.

Fabricating a prosthesis that is optimally tailored to the patient demands technical knowledge of orthopaedic technology. The alignment and adjustments may only be carried out by a prosthetist. The alignment and adjustments are mainly completed in the following steps:

1. Bench alignment
2. Static alignment
3. Dynamic trial fitting

These steps are first carried out with a test prosthesis in order to determine the best possible combination and positioning of the prosthetic components relative to each other. Once all settings have been tailored to the patient, the definitive prosthesis is produced. The same steps are repeated for this purpose.

In dynamic trial fitting, the patient must be given sufficient time to become familiar with the functions of the prosthesis, thereby learning safe use for everyday life.

5.2 Information on Using a Cosmetic Cover

CAUTION

Use of talcum

Risk of injury, damage to the product due to lack of lubrication

- Do not use talcum on the product or other prosthesis components.

- Only use the foam cover approved by the manufacturer (see section "Combination Possibilities" – see page 20).
- To reduce friction and to eliminate noise, apply 519L5 Silicone Spray directly onto the contact surfaces of the cosmetic foam cover.
- After fabricating the cosmetic cover, repeat the dynamic trial fitting since the cover affects the adjustments of the prosthesis (e.g. damping).

5.3 Bench Alignment

Use 2Z11=KIT

INFORMATION: The connection area of the prosthetic knee joint can be protected against scratches with the protective film of the 2Z11=KIT during alignment in the workshop and testing in the trial fitting area.

- ▶ Use the protective film as illustrated in the accompanying document for the 2Z11=KIT.
- ▶ Remove the protective film before the patient leaves the fitting area.

Bench alignment process	
+ = anterior placement / - = posterior placement (in relation to the alignment reference line)	
Item	Grundaufbau
	<p>Required equipment: Alignment apparatus (e.g., 743L200 L.A.S.A.R. Assembly or 743A200 PROS.A. Assembly), 743A80 50:50 Gauge, 743S12 Heel Height Measuring Apparatus</p> <p>Position the prosthetic foot in the alignment apparatus in accordance with the following values:</p>
①	<p>Heel height: effective heel height (x) + 5 mm</p>
②	<p>a-p positioning, middle of the prosthetic foot to the sagittal alignment reference line ①: +30 mm</p> <p>Exterior foot rotation: approx. 7°</p> <p>m-l positioning, middle of the prosthetic foot to the frontal alignment reference line ⑤: 0 mm medial</p>
③	<p>Check the prosthetic knee joint according to the section "Checking the Factory Settings" (see page 29). Position the alignment reference point ② (upper, front rotation axis) of the prosthetic knee joint in the alignment apparatus according to the following values: Height of the alignment reference point in the alignment apparatus: MTP-floor measurement + 20 mm</p> <p>a-p positioning of the alignment reference point to the alignment reference line: 0 mm</p> <p>m-l positioning, centre of the prosthetic knee joint to the frontal alignment reference line ⑤: 0 mm</p> <p>External rotation of the prosthetic knee joint: approx. 5°</p>
④	<p>INFORMATION: If the length-adjustable 2R45=S Tube Adapter (max. body weight 100 kg) is used in the alignment of the test prosthesis, the required length for the definitive prosthesis is easier to determine. When using the included tube adapter, note the section "Adapting the Length of the Tube Adapter" (see page 25). Connect the prosthetic knee joint to the prosthetic foot and chosen adapters according to their instructions for use and the section "Installing the Tube Adapter" (see page 25).</p>
⑤	<p>Position the 743A80 50:50 Gauge on the lateral side of the prosthetic socket and mark a point respectively in the distal and proximal area. Using these points, mark the lateral centre line of the prosthetic socket. Mark the tuberosity height as the tuberosity reference point ④ on the centre line. Clearly mark the lateral socket reference point ③ 30 mm proximal to the tuberosity reference point. Position the tuberosity reference point ④ of the prosthetic socket at the height of the tuberosity-floor measurement.</p>

Bench alignment process	
+ = anterior placement / - = posterior placement (in relation to the alignment reference line)	
Item	Grundaufbau
6	<p>CAUTION! The correct position and flexion of the prosthetic socket is crucial for the safety of the patient and functionality of the prosthetic knee joint, and avoids premature wear.</p> <p>INFORMATION: The 4R118 Adapter Plate (max. body weight 125 kg) can be used for further posterior placement.</p> <p>Position the prosthetic socket so the lateral socket reference point ③ lies on the sagittal alignment reference line ①.</p> <p>Adjust the socket flexion angle β depending on the hip flexion contracture (angle α) according to the following specifications:</p> <p>If the angle $\alpha = 0^\circ$, the optimum value of the angle $\beta = 3 \text{ to } 5^\circ$ to the alignment reference line.</p> <p>If the angle $\alpha > 0^\circ$, the optimum value of the angle $\beta = \alpha + 5 \text{ to } 10^\circ$.</p>
7	<p>To fabricate the prosthetic socket for the 3R106-Pro=KD, observe the section "Fabricating the prosthetic socket" (see page 26).</p> <p>Position the lamination anchor or socket adapter so that all prosthetic components between the socket and prosthetic knee joint can be connected to each other correctly.</p> <p>Join the prosthetic socket and the prosthetic knee joint using the adapters you selected.</p> <p>Be sure to follow the adapter instructions for use during adjustment and installation.</p> <hr/> <p>Position the 743A80 50:50 Gauge on the anterior side of the prosthetic socket and mark a point respectively in the distal and proximal area.</p> <p>Using these points, mark the anterior centre line of the prosthetic socket.</p> <p>Clearly mark the anterior socket reference point ⑥ (intersection of the socket edge and anterior centre line).</p> <p>Position the prosthetic socket so the anterior socket reference point ⑥ lies on the alignment reference line ⑤.</p> <p>Adjust the socket flexion angle γ according to the patient's adduction angle.</p>

5.3.1 Shortening the Tube Adapter

CAUTION

Incorrect processing of tube
 Fall due to damage to the tube

- ▶ Do not clamp the tube into a vice.
- ▶ To shorten the tube, use only a tube cutter or a cutting device.

> **Tools required:**

719R3 Tube Cutter, 718R1 Tube Deburrer

- 1) **CAUTION! The length of the tube must not be less than the required length for the patient, since length compensation in the insertion area of the prosthetic knee joint is prohibited.**

Use the tube cutter to shorten the tube to the required length (Rohr schneiden).

- 2) Carefully deburr the inside and outside of the tube with the tube deburrer (Rohr entgraten).
- 3) Dispose of the remainder.

5.3.2 Installing the Tube Adapter

CAUTION

Incorrect mounting of the tube
 Risk of injury due to breakage of load-bearing components

- ▶ Slide the tube all the way to the stop in the intended prosthetic component when mounting.

CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

Installation on the Prosthetic Knee Joint

> **Required tools and materials:**

710D4 Torque Wrench, degreasing cleaner (verify material compatibility)

- 1) Loosen the cap screw by **approx. 2 turns** (Rohradapter festschrauben).
- 2) Clean the contact surfaces of the prosthetic knee joint and tube adapter using a degreasing cleaning agent.
- 3) Slightly rotate the tube as you insert it into the prosthetic knee joint to the stop (Rohradapter einsetzen).
- 4) Using the torque wrench, tighten the cap screw (**10 Nm**) (Rohradapter festschrauben).

Installation on the Prosthetic Foot

> **Required tools and materials:**

710D4 Torque Wrench, 636K13 Loctite®

- 1) Position the pyramid receiver of the tube adapter on the pyramid adapter of the prosthetic foot.
- 2) **For final assembly:**
Use Loctite® to secure the set screws of the pyramid receiver.
- 3) Screw in the set screws.
- 4) Use the torque wrench to tighten the set screws (**10 Nm**).
- 5) Replace any set screws that are protruding or recessed too much with suitable ones (see selection table).
- 6) **For final assembly:**
Use the torque wrench to retighten the set screws (**15 Nm**).

Selection table for set screws

Reference number	Length (mm)
506G3=M8X12-V	12
506G3=M8X14	14
506G3=M8X16	16

5.3.3 Fabricating the prosthetic socket

INFORMATION

The instructions for fabricating the prosthetic socket are intended only for the 3R106-Pro=KD Prosthetic Knee Joint.

5.3.3.1 Fabricating the Test Socket

Laminating the Prosthetic Socket

> **Required materials:**

640Z5* Plaster Insulating Cream, 623T3* Perlon Stockinette, 616G13* Woven Fibreglass Stockinette, 616G12* Carbon Fibre Cloth, 617H19* Orthocryl Lamination Resin, 617H17* Orthocryl soft, 4G70 Lamination Anchor

- 1) Apply plaster insulating cream to the plaster positive.
- 2) Cut a piece of Perlon stockinette to twice the length of the plaster positive and pull on one half.

- 3) Twist the Perlon stockinette on the distal end and pull on the other half.
- 4) Cut a piece of woven fibreglass stockinette to the corresponding length and pull on the length equal to two-thirds of the prosthetic socket.
- 5) Tie off the woven fibreglass stockinette on the distal end and pull on the length equal to half the prosthetic socket.
- 6) **INFORMATION: For stability, a different alignment of the carbon fibres and an overlap of 3 cm to the edge of the lamination anchor are important.**
Cut two pieces of carbon fibre cloth to size as a mat according to the edge of the lamination anchor.
- 7) Position the two pieces on the distal area of the plaster positive.
- 8) Cut a piece of Perlon stockinette to twice the length of the plaster positive and pull on one half.
- 9) Twist the Perlon stockinette on the distal end and pull on the other half.
- 10) **INFORMATION: The double casting process is used for lamination.**
Cast the distal section of the prosthetic socket ($\frac{2}{3}$ of the length) with Orthocryl lamination resin and allow it to cure.
- 11) Cut a piece of Perlon stockinette to twice the length of the plaster positive and pull on one half.
- 12) Twist the Perlon stockinette on the distal end and pull on the other half.
- 13) Cast the prosthetic socket with Orthocryl soft and allow it to cure.

Applying the Lamination Anchor

> Required materials:

636K8=20x2x10 Plastaband, Adhesive Putty (consisting of: 639A1=1 Talcum, 617H21* Orthocryl Sealing Resin, 617P37=0.150 Hardener Powder), 627B2* Monofilament Adhesive Tape

- 1) Bend the lamination anchor to fit the prosthetic socket.
- 2) Fill the coupling opening of the lamination anchor with Plastaband.
- 3) Position the lamination anchor on the prosthetic socket and glue it in place with adhesive putty.
- 4) Secure the glued lamination anchor to the prosthetic socket with adhesive tape.
- 5) Bring the prosthetic knee joint to full flexion. If necessary, form the required contact surface on the prosthetic socket with putty. It is also possible to glue a stop bumper of Pedilin onto the outer laminate.

5.3.3.2 Fabricating the definitive socket

> Required materials:

4x501T1=M5X16 Slotted Truss Head Screw, 4x 502R1=M5X16 Two Hole Nut, 616G12* Carbon Fibre Woven Fabric, 623T3* Perlon Stockinette, 616G13* Woven Fibreglass Stockinette, 617H19* Orthocryl Lamination Resin, 617H17* Orthocryl soft, 4G70 Lamination Anchor

- 1) Secure the connection between the lamination anchor and prosthetic socket with the help of the truss head screws and two hole nuts.
- 2) **INFORMATION: For stability, a different alignment of the carbon fibres and an overlap of 3 cm to the edge of the lamination anchor are important.**
Cut two pieces of carbon fibre cloth to size as a mat according to the edge of the lamination anchor.
- 3) Position the two pieces on the lamination anchor.
- 4) Cut a piece of Perlon stockinette to twice the length of the plaster positive and pull on one half.
- 5) Twist the Perlon stockinette on the distal end and pull on the other half.
- 6) Cut a piece of woven fibreglass stockinette to the corresponding length and pull on the length equal to two-thirds of the prosthetic socket.
- 7) Tie off the woven fibreglass stockinette on the distal end and pull on the length equal to half the prosthetic socket.

- 8) **INFORMATION: The double casting process is used for lamination.**
Cast the distal section of the prosthetic socket (2/3 of the length) with Orthocryl lamination resin and allow it to cure.
- 9) Cut a piece of Perlon stockinette to twice the length of the plaster positive and pull on one half.
- 10) Twist the Perlon stockinette on the distal end and pull on the other half.
- 11) Cast the prosthetic socket with Orthocryl soft and allow it to cure.

5.4 Static Alignment

Static alignment process	
+ = anterior placement / - = posterior placement (in relation to the load line)	
Item	Statischer Aufbau
	Required equipment: 743L100 L.A.S.A.R. Posture
①	To determine the sagittal load line ①, position the patient on the L.A.S.A.R. Posture as follows: <ul style="list-style-type: none"> • Prosthetic foot (with shoe) on force measurement plate (sufficient load: > 35% body weight) • Other foot (with shoe) on height compensation plate • The tips of the shoes are in one line
②	Optimise the static alignment solely by changing the plantar flexion. Adjustment is performed only with the anterior and posterior set screws of the tube adapter on the prosthetic foot.
	a-p positioning of the alignment reference point ② (upper, anterior rotation axis) to the sagittal load line ①: -35 mm
	m-l positioning, middle of the prosthetic foot to the frontal load line ④: 1 to 2 cm medial
	m-l positioning, centre of the prosthetic knee joint to the frontal load line ④: 1 to 2 cm medial
	m-l positioning of the prosthetic socket: Positioning of the anterior superior iliac spine (Spina iliaca anterior superior) ③ to the frontal load line ④: 1 to 2 cm lateral

5.5 Dynamic Trial Fitting

⚠ CAUTION

Adjusting the Settings
Fall due to incorrect or unfamiliar settings

- ▶ Only adapt the settings to the patient gradually.
- ▶ Explain the effects of the adjustments on the use of the prosthesis to the patient.

During dynamic trial fitting, the alignment and settings of the prosthesis are checked and adapted for optimum walking according to the needs and abilities of the patient.

The patient has to learn the safe use of the prosthesis through intensive training.

The following subsections describe the adjustment possibilities for adapting the product to the patient.

The list below provides an overview of the sequence for working through the subsections:

- Checking the Factory Settings
- Setting Flexion Damping
- Setting Extension Damping

5.5.1 Checking the Factory Settings

Function	Means of adjustment	Factory setting	Meaning
Flexion damping	Adjustment valve "F" (see fig. 4 - ①)	From the left-hand stop, turned to the right (number of turns: 2.0)	Valve slightly closed – low damping as the starting point for adjustment to the patient
Extension damping	Adjustment valve "E" (see fig. 4 - ②)	Turned all the way to the left	Valve fully open – minimal damping as the starting point for adaptation to the patient (safety setting)

> **Tools required:**

710H10=2x3 Adjustment Wrench

- 1) Check whether the settings match the factory settings (see table).
- 2) In case of deviations, restore the factory settings using the adjustment wrench.

5.5.2 Setting Flexion Damping

Function	Means of adjustment	Setting	Meaning
Flexion damping	Adjustment valve "F" (see fig. 4 - ①)	Rotation to the right	Increase damping
		Rotation to the left	Reduce damping

> **Tools required:**

710H10=2x3 Adjustment Wrench

▶ **CAUTION! The prosthetic foot must not swing through too far.**

Adjust flexion damping using the adjustment wrench so that the prosthetic foot swings through sufficiently according to the needs of the patient.

5.5.3 Setting Extension Damping

Function	Means of adjustment	Setting	Meaning
Extension damping	Adjustment valve "E" (see fig. 4 - ②)	Rotation to the right	Increase damping
		Rotation to the left	Reduce damping

> **Tools required:**

710H10=2x3 Adjustment Wrench

▶ **CAUTION! The prosthetic knee joint has to reach full extension, even at slow walking speeds. If the 513D83=1.1X12.9X66 Compression Spring, weak is installed and the extension assist effect is not sufficient for the patient, the strong compression spring has to be installed again (section "Installing the compression spring" –).**

Using the adjustment wrench, adjust extension damping so that the prosthetic knee joint does not swing too hard against the extension stop.

5.5.4 Optional: Replacing the Compression Spring

If required, the installed strong compression spring can be replaced with the weak compression spring or omitted entirely (see the section "Scope of Delivery" – see page 22). **To prevent injuries due to falling**, this is only permitted if the patient is able to bring the prosthetic knee joint to full extension prior to heel strike in all situations.

To prevent damage to the product, the centring disc ② must be removed from the prosthetic knee joint using suitable pliers (see fig. 6) if the compression spring is omitted entirely.

> **Tools required:**

Flat screwdriver

- 1) Keep the prosthetic knee joint vertical during assembly and disassembly (see fig. 6).
- 2) Insert the flat screwdriver ④ into the groove ③ of the bushing ① (see fig. 5 – ①).
- 3) With the flat screwdriver, first push towards the centre until the cam of the bushing no longer engages the groove of the clamp bracket ② (see fig. 5 – ②).
- 4) Slide the bushing out of the clamp bracket using the flat screwdriver (see fig. 5 – ③).

- 5) Remove the compression spring ③ from the bushing (see fig. 6).
- 6) Insert the other compression spring in the bushing (see fig. 6).
- 7) If the centring disc ② fell out during disassembly, set the centring disc onto the compression spring (see fig. 6).
- 8) Bring the prosthetic knee joint to 90° flexion.
- 9) Slide the bushing as far as possible into the clamp bracket.
- 10) **INFORMATION: The edge of the bushing must be fully in contact with the clamp bracket.**
Press together the two ends of the bushing and slide the bushing into the clamp bracket to the stop.
- 11) **INFORMATION: The tube adapter will only slide into the clamp bracket if the cam of the bushing has fully engaged in the groove of the clamp bracket. The contact surface of the bushing to the tube adapter must be free of grease.**
Use a finger to push the cam of the bushing into the groove in the clamp bracket.
- 12) Check the settings of the prosthetic knee joint and reaching extension with a walking test.

5.6 Finishing the Prosthesis

> Required tools:

Torque wrench

- 1) Tighten all screws of the prosthetic components to the specified installation torque values.
- 2) Verify proper functioning of the prosthesis.

6 Use

6.1 Information for use

NOTICE

Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

INFORMATION

At a flexion angle > 110°, the compression spring automatically moves the prosthetic knee joint against the flexion stop. This function is activated in particular when kneeling.

- ▶ **Inform the patient.**

6.2 Cleaning

⚠ CAUTION

Use of unsuitable cleaning agents or disinfectants

Impairment of functionality and damage due to incorrect cleaning agents or disinfectants

- ▶ Only clean the product with the approved cleaning agents.
- ▶ Only disinfect the product with the approved disinfectants.
- ▶ Observe the instructions for cleaning and care.

INFORMATION

- ▶ Clean the product if it gets dirty.

- ▶ For disinfection, only use disinfectants that do not affect the materials of the product. For further information, please consult the manufacturer.
- ▶ Observe the cleaning instructions for all prosthetic components.

Cleaning light soiling

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

Cleaning heavier soiling

> **Required materials:** Cleaning cloth, 634A58 Isopropyl Alcohol, soft cloth

- 1) **NOTICE! Verify material compatibility! Do not remove lubricant from bearing points!**
Clean the product with a cleaning cloth and isopropyl alcohol.
- 2) Dry the product with a cloth.

Cleaning with Disinfectants

> **Required materials:** colourless, alcohol-free disinfectant (verify material compatibility!), soft cloth

- 1) Disinfect the product with the disinfectant.
- 2) Dry the product with a cloth.
- 3) Allow to air dry in order to remove residual moisture.

7 Maintenance

⚠ CAUTION

Failure to follow the maintenance instructions

Risk of injuries due to changes in or loss of functionality and damage to the product

- ▶ Observe the following maintenance instructions.

- ▶ **NOTICE! Do not lubricate and grease the prosthetic joint.**
- ▶ **NOTICE! Repair work must be performed exclusively by manufacturer service.**
- ▶ Arrange regular maintenance intervals with the patient depending on the level of use.
- ▶ Following an individual period for the patient to get accustomed to the prosthesis, check the settings of the prosthetic joint and adapt them to the patient's requirements again as needed.
- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.
- ▶ As part of the safety inspections, inspect the prosthetic joint for wear and proper functionality. Special attention should be paid to movement resistance, bearings and abnormal noises. Full flexion and extension must be ensured at all times. Readjust as required.

8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Legal information

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

10 Technical data

Reference number	3R106-Pro	3R106-Pro=KD	3R106-Pro=ST
Weight [g]	885	910	915
System height [mm]	163	187	181
Proximal system height up to the alignment reference point [mm]	-7	16	11
Distal system height up to the alignment reference point [mm]	170		
Knee flexion angle	175°		
Proximal connection	Pyramid Adapter	Lamination Anchor	Threaded Connector (M26x1.5)
Distal connection	Tube Clamp (Ø 30 mm)		
Max. body weight [kg]	125		
Mobility grade	2 + 3		

Reference number	2R49
Weight [g]	240
Material	Aluminium
Diameter [mm]	30
Min. system height [mm]	97
Max. system height [mm]	432
Distal connection	Pyramid Receiver
Max. body weight [kg]	125

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