
17B44

DE	Gebrauchsanweisung	3
EN	Instructions for use	9
FR	Instructions d'utilisation	15
IT	Istruzioni per l'uso	22
ES	Instrucciones de uso	28
PT	Manual de utilização	35
NL	Gebruiksaanwijzing	41
SV	Bruksanvisning	48
DA	Brugsanvisning	53
NO	Bruksanvisning	59
PL	Instrukcja użytkowania	65
HU	Használati utasítás	71
CS	Návod k použití	77
TR	Kullanma talimatı	83
EL	Οδηγίες χρήσης	89
RU	Руководство по применению	96

1

INFORMATION

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INFORMATION

Date of last update: 2020-04-02

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the processing of the 17B44 orthotic knee joint with Swiss locking mechanism.

2 Product description

2.1 Components/design

Included in the scope of delivery (see fig. 1)		
Item:	Component	Reference number
1	Lock lever	17Y57
2	Lock screw	501A23
3	Bolt	
4	Lock plate	
5	Joint screw	501S55=M6x20x16

Included in the scope of delivery (see fig. 1)		
Item:	Component	Reference number
6	Countersunk oval-head screw	501T7=7.5x9xM5
7	Lower joint section	
8	Knee cap ring plate	17Y56
9	Upper joint section	

Not included in scope of delivery (see fig. 1)			
No.	Qty.	Description	Reference number
10	2	Side bar	17B38=*, 17B39=*, 17B40=*, 17B41=*, 17B84=*
not illustrated	2	Lamination bar	17Y128=*, 17Y129=*

2.2 Design

The locking mechanism consists of a centred lock lever and a toothed lock plate. A toothed bolt is installed in the lower joint section. The orthotic joint angle can be adjusted in 8° increments by repositioning the lock plate on the upper joint section.

3 Intended use

3.1 Indications for use

The orthotic joint is intended **exclusively** for use as a pair for orthotic fittings of the lower limbs on one patient.

INFORMATION

Increased wear and tear

The orthotic knee joints are intended only for use in combination with free-motion or foot lifting orthotic ankle joints. Wear occurs more rapidly in combinations involving orthotic ankle joints with a dorsal stop.

INFORMATION

With flexion adjustment > 16°, the knee joint is to be used only for positioning orthoses.

3.2 Area of application

Approved for a body weight of **up to 100 kg**.

3.3 Indications

Paresis or total paralysis of the leg muscles with existing extension deficit in the knee joint.

Indications must be determined by the physician.

4 Safety

4.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

4.2 Safety instructions

 **CAUTION**

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

 **CAUTION**

Incorrect alignment or assembly

Injuries due to changes in or loss of functionality

- ▶ The product may only be installed by trained, qualified personnel.
- ▶ Observe the alignment and assembly instructions.

 **CAUTION**

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

 **CAUTION**

Improper handling

Risk of falling due to accidentally opening the lock

- ▶ Prior to handing over the orthosis, ensure that the locks are engaged securely and in parallel.
- ▶ Avoid accidentally opening the locks.
- ▶ Instruct the patient in the correct usage of the system joints.

⚠ CAUTION

Clothing getting caught in the orthotic joint

Risk of injury due to loss of orthotic joint functionality

- ▶ Make sure your clothing does not get caught in the orthotic joint.

⚠ CAUTION

Excessive strain due to use on more than one patient

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

NOTICE

Thermal overloading of the product

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at temperatures above **300 °C (570 °F)**.
- ▶ Prior to thermal treatment, remove all temperature-critical components (such as plastic parts).

NOTICE

Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below **-10 °C (14 °F)** or above **+60 °C (140 °F)** (e.g. sauna, excessive sunlight, drying on a radiator).

5 Preparation for use

INFORMATION

Corrosion protection

The side bars and orthotic joints are resistant to corrosion. Fine grinding and polishing increase corrosion protection. Use sintering powder (e. g. 618T40) for coating.

Adjusting the orthotic joint angle

- 1) Unlock the orthotic joint with the lock lever.
- 2) Unscrew the joint screw.
- 3) Loosen the lock plate from the bolt.
- 4) Adjust the orthotic joint angle by turning the lock plate.
- 5) Insert the joint screw into the lock plate and tighten it (tightening torque: **12 Nm**).

Gluing in bars

Required materials: 636W18 special glue for bars, 636W19 hardener, optional: 636W28 special adhesive, degreasing cleaner

The bars must be glued in the insertion zones to ensure stability.

- 1) Clean the connecting surfaces with a degreasing cleaner.
- 2) Mix the special adhesive and hardener.
- 3) Apply the mixture to the insertion zones.
- 4) Insert the bars.
- 5) Secure with the enclosed screws.
- 6) Allow to cure for at least **4 hours**.

INFORMATION: Final bonding strength will be reached after 16 hours.

Interconnecting the lock levers

For simultaneous operation of the two 17Y57 lock levers, the lock levers can be interconnected.

- 1) Weld or solder the 17Y53 connection tube to the lock lever.
- 2) Verify that the two locks release and engage simultaneously.
- 3) Connect the lock levers with a traction element (e. g. elastic strap) to enable engaging.

6 Cleaning

Clean the orthosis joint promptly after contact with water containing salt, chlorine or soap, or if it gets dirty.

- 1) Rinse the orthosis joint with clean, fresh water.

- 2) Dry the orthosis joint with a cloth or allow it to air-dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

7 Maintenance

⚠ CAUTION

Use of an incorrect lock

Injuries from falling due to an incorrect lock being installed

- ▶ Only use locks of the same size when reinstalling or replacing.

INFORMATION

The product may be exposed to increased loads by the patient.

- ▶ Shorten the maintenance intervals according to the expected loads.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

As part of maintenance, inspect the orthotic joints for wear and proper functionality.

Spare parts are listed in the section "Components/design".

In case of wear and tear on the bolt and lock plate, send the orthotic knee joint to the Ottobock Service Department for repair.

Only use 633F7 special lubricant.

Replacing the lock lever and lock screw

- 1) Loosen the lock screw with a screwdriver and remove it.
- 2) Remove the lock lever.
- 3) Insert the replaced lock lever into the upper joint section.

INFORMATION: When replacing the lock lever, replace the lock screw as well.

- 4) Insert the lock screw and secure with Loctite 241.
- 5) Following installation, adjust the lock lever according to the needs of the patient.

INFORMATION: Make adjustments only on the lock lever. Do not modify the lock plate and bolt.

8 Disposal

Dispose of the product in accordance with national regulations.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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