# ottobock.

## 8K26, 8K27

<b>DE</b>	Gebrauchsanweisung .....	3
<b>EN</b>	Instructions for use .....	8
<b>FR</b>	Instructions d'utilisation .....	13
<b>IT</b>	Istruzioni per l'uso .....	18
<b>ES</b>	Instrucciones de uso .....	23
<b>PT</b>	Manual de utilização.....	29
<b>NL</b>	Gebruiksaanwijzing.....	34
<b>SV</b>	Bruksanvisning.....	40
<b>DA</b>	Brugsanvisning.....	44
<b>NO</b>	Bruksanvisning.....	49
<b>PL</b>	Instrukcja użytkowania.....	54
<b>HU</b>	Használati utasítás .....	59
<b>CS</b>	Návod k použití.....	64
<b>RO</b>	Instrucțiuni de utilizare.....	68
<b>HR</b>	Upute za uporabu .....	74
<b>TR</b>	Kullanma talimatı.....	79
<b>EL</b>	Οδηγίες χρήσης .....	84
<b>RU</b>	Руководство по применению .....	89
<b>JA</b>	取扱説明書.....	95
<b>ZH</b>	使用说明书.....	99

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## Explanation of symbols

English

**CAUTION** Warnings regarding possible risks of accident or injury.

**NOTICE** Warnings regarding possible technical damage.

### INFORMATION

Last update: 2015-12-02

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

The scope of delivery is shown on the front page.

## 1 Single components

### 1.1 Single components

For available single components refer to the catalogue.

### 1.2 Accessories

Select the glove according to the hand size.

8S4=*	Prosthetic Glove
8S5=*	Prosthetic Glove
8S6=*	Prosthetic Glove
8S4N=*	MovoSkin Natural
8S5N=*	MovoSkin Natural
8S6N=*	MovoSkin Natural
640F12=	Special Cleaner
640F13=	Pump Sprayer for 640F12

## **Accessories for long below-elbow and through-wrist residual limb**

9S96=40	Chassis	(for hand size L / R 6 ¾)
9S184=*	Chassis	(for hand sizes L / R 7 ¼, L / R 7 ¾, and L / R 8)
506G4=M3x5	Set Screw	(for hand size L / R 6 ¾)
506G4=M4x5	Set Screw	(for hand sizes L / R 7 ¼, L / R 7 ¾, and L / R 8)
9E94=*	Lamination Ring	(for hand sizes L / R 6 ¾, L / R 7 ¼, L / R 7 ¾ and L / R 8)

The accessories are not included in the delivery. They must be ordered separately.

## **2 Description**

### **2.1 Indications for use**

The Ottobock System Hands –voluntary closing– are to be used solely for the exoprosthetic fitting of the upper extremities in combination with Ottobock arm components.

### **2.2 Field of application**

The System Hands –voluntary closing– can be used for any residual limb length. In individual cases it must be tested whether the wearer of the arm prosthesis is able to completely close the hand with a harness control system.

### **2.3 Safety instructions**

Please make your patients aware of the following safety instructions.

#### **⚠ CAUTION**

**Risk of injury as a result of excessive strain.** The Ottobock prosthetic components described in these instructions for use have been developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's safety!

Should prosthetic components be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the components inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Service.

### CAUTION

**Risk of accident when driving a motor vehicle.** The ability of prosthesis users to drive a vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities.

All persons are required to observe their country's national and state driving laws when operating motor vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorized test centre.

For maximum safety and convenience, Ottobock recommends that, at the very least, a specialist evaluate the need for any adaptations to the car (such as by installing a steering fork, automatic shift). Risk-free driving should be ensured even when the prosthesis is not functioning.

### NOTICE

**Risk of corrosion.** Prosthetic components must not be exposed to corrosive elements such as fresh water, salt water, and acids.

Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

### NOTICE

**Damage caused by inappropriate environmental conditions.** Prosthetic components must not be subjected to intense smoke, dust, vibrations, shocks, or high temperatures. Do not allow debris or liquids to get into the prosthetic components.

Failure to follow these instructions can lead to malfunctioning and damage to the prosthesis.

## NOTICE

**Damage caused by improper cleaning agents.** Do not use any aggressive cleaning agents since they could cause damage to bearings, seals, and plastic parts.

### 2.4 Function

The System Hands –voluntary closing– differ from each other, among other things, by the inner pull and outer pull. The inner pull is provided with the cable control through the palm of the hand. The outer pull is provided with the cable control through the back of the hand.

The Ottobock System Hands –voluntary closing– are available in the following designs:

8K26=L/R\* with metric thread M12x1.5 with outer pull

8K27=L/R\* with inch thread ½"-20 with outer pull

For the various sizes see table in section 3.

The System Hands –voluntary closing– belong to the cable-activated prostheses and are closed by a harness control system (active). The grip force is increased by intensifying the pull. Further pull releases the lock and spring tension opens the hand.

The System Hands –voluntary closing– perform the three-finger grip with the index, middle finger, and thumb.

To get a harmonious picture of the closing and opening movement of the entire hand, the ring finger and small finger are coupled to the middle finger via the inner hand.

The grip force needed for opening can be changed by adjusting the spiral spring. For this purpose, the System Hand –voluntary closing– is to be sent to an Ottobock Service centre.

### 3 Technical data

Size	6 ¾	7 ¼	7 ¾	8
Weight [g]	250	340	370	380
Total height [mm]	125	152	155	157
Thumb length [mm]	105	132	135	137

## **4 Handling**

### **4.1 Assembly**

The System Hands –voluntary closing– are to be mounted to an Ottobock wrist unit.

For long below-elbow and through-wrist residual limbs, the chassis with threaded stud must be replaced by the chassis with 9E94=\* Lamination Ring as follows:

9S96=40 for hand size L/R 6 ¾

9S184=\* for hand sizes L/R 7 ¼, L/R 7 ¾, and L/R 8

The parts are attached with the 506G4=\* Set Screw (Fig. 1).

The 640F18 Donning Spray facilitates the application of the inner hand on the System Hand.

To avoid any damage, use the 711M1 Hand Holding Tool for holding the System Hand while pulling on the glove. The 711M2 Mounting Plate is used for attaching the Hand Holding Tool to a workbench (Fig. 2).

### **4.2 Combination possibilities**

The Ottobock System Hands –voluntary closing– can be combined with the adapters of the field of modular arm prostheses and with the Ottobock wrist units.

It may be necessary to select a suitable adapter from the wrist unit accessories.

## **5 Legal information**

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### **5.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### **5.2 CE conformity**

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device

according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH