
# ottobock.

## 10V8, 10V30

|           |                                  |     |
|-----------|----------------------------------|-----|
| <b>DE</b> | Gebrauchsanweisung .....         | 3   |
| <b>EN</b> | Instructions for use .....       | 9   |
| <b>FR</b> | Instructions d'utilisation ..... | 14  |
| <b>IT</b> | Istruzioni per l'uso .....       | 21  |
| <b>ES</b> | Instrucciones de uso .....       | 27  |
| <b>PT</b> | Manual de utilização .....       | 33  |
| <b>NL</b> | Gebruiksaanwijzing .....         | 39  |
| <b>SV</b> | Bruksanvisning .....             | 46  |
| <b>DA</b> | Bruksanvisning .....             | 51  |
| <b>NO</b> | Bruksanvisning .....             | 57  |
| <b>FI</b> | Käyttöohje .....                 | 62  |
| <b>PL</b> | Instrukcja użytkowania .....     | 68  |
| <b>HU</b> | Használati utasítás .....        | 74  |
| <b>CS</b> | Návod k použití .....            | 80  |
| <b>RO</b> | Instrucțiuni de utilizare .....  | 86  |
| <b>TR</b> | Kullanma talimatı .....          | 93  |
| <b>EL</b> | Οδηγίες χρήσης .....             | 99  |
| <b>RU</b> | Руководство по применению .....  | 105 |
| <b>JA</b> | 取扱説明書 .....                      | 111 |
| <b>ZH</b> | 使用说明书 .....                      | 117 |

## INFORMATION

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## INFORMATION

Last update: 2021-02-15

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

## Explanation of Symbols

English

**⚠ WARNING** Warnings regarding possible risks of severe accident or injury.

**⚠ CAUTION** Warnings regarding possible risks of accident or injury.

**NOTE** Warnings regarding possible technical damage.

**INFORMATION** Additional information on the fitting/use.

## 1 Single components

### 1.1 Single components

For available single components refer to the catalogue.

### 1.2 Accessories

*10V8*

10A43=\* Attachment Plate with stud

10V15 Connector

10V26 Flexion Device for later installation

*10V30*

10A56=\* Attachment Plate with short stud

The accessories are not included. They must be ordered separately.

## 2 Description

### 2.1 Indications for use

The 10V8/10V30 Ottobock Wrist Units – ratchet type rotation – are to be used **solely** for the exoprosthetic fitting of the upper extremities in combination with Ottobock arm components.

### 2.2 Field of application

The 10V8/10V30 Ottobock Wrist Units – ratchet type rotation – can be used in combination with an Ottobock Cable-Activated Hook, an All-Purpose Hook or an Ottobock System Hand.

The 10V8 is the ratchet type rotation wrist unit for standard fittings.

The 10V30 is used for long forearm residual limbs, in cases where a fitting with the 10V8 would result in an excessive arm length.

### 2.3 Lifetime

Expected lifetime: 5 years

The design, manufacturing and requirements for the intended use of the product are based on the expected lifetime.

### 2.4 Safety instructions

Please make your patients aware of the following safety instructions:

#### **WARNING**

**Risk of accident when driving a motor vehicle.** The ability of prosthesis users to drive a vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities.

All persons are required to observe their country's national and state driving laws when operating motor vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorized test centre.

For maximum safety and convenience, Ottobock recommend that, at the very least, a specialist evaluate the need for any adaptations to the car (such as by installing a steering fork, automatic shift). Risk-free driving should be ensured even when the prosthesis is not functioning.

**⚠ CAUTION**

**Risk of injury as a result of improper manipulation.** Inadvertent activation of the release lever can make the terminal device come loose.

**⚠ CAUTION**

**Risk of injury as a result of excessive strain.** The Ottobock prosthetic components described in these instructions for use have been developed for activities of daily living and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's safety!

Should prosthetic components be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the components inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Service Centre.

**NOTICE**

**Risk of corrosion.** Prosthetic components must not be exposed to corrosive elements such as fresh water, salt water, and acids.

Using this medical product in such environmental conditions will render all claims against Ottobock null and void.

**NOTICE**

**Damage caused by inappropriate environmental conditions.** Prosthetic components must not be subjected to intense smoke, dust, vibrations, shocks, or high temperatures. Do not allow debris or liquids to get into the prosthetic components. Failure to follow these instructions can lead to malfunctioning and damage to the prosthesis.

**NOTICE**

**Damage caused by improper cleaning agents.** Do not use any aggressive cleaning agents, since they could cause damage to bearings, seals and plastic parts.

## 2.5 Function

The Ottobock Wrist Units – ratchet type rotation – are available in the following designs:

- 10V8 Wrist Unit – ratchet type rotation
- 10V30 Wrist Unit – ratchet type rotation, short

The ratchet type rotation wrist units feature two levers.

The locking lever A locks rotation of the Ottobock System Hand or Ottobock Cable-Activated Hook (Fig. 1, lever A). Actuating the lever once again will unlock rotation. Now the terminal device can be rotated to one of the 12 different ratchet positions.

The attachment plate automatically locks in place, when the plate's stud audibly engages with the terminal stop when pushing in.

Actuating the release lever B (Fig. 1, lever B) releases the attachment plate with stud from the wrist unit.

Connecting the 10V26 Flexion Device to the 10V8 Wrist Unit – ratchet type rotation – results in two other movement planes.

When using the 10V26 Flexion Device, the flexion of the 10V8 Wrist Unit – ratchet type rotation – can be locked with the lever (Fig. 2).

## 3 Technical data

| Article number    | 10V8 | 10V30 |
|-------------------|------|-------|
| Weight [g]        | 120  | 100   |
| Total height [mm] | 25   | 19    |

## 4 Handling

### 4.1 Assembly

Clean the threads of the 10A43 and 10A56 Attachment Plates with stud and of the terminal device with acetone and let air dry. Next, apply a little bit of 636K14 Loctite 601 to the thread of the terminal device. Now screw the stud onto the thread of the terminal device up to the stop.

## NOTICE

**Risk of damage as a result of improper handling of the glue.** To ensure a solid bond, the respective areas must be degreased with acetone before applying the adhesive.

Once you have degreased the areas, do not touch them anymore and immediately apply the adhesive to the surfaces after the acetone has dried.

## INFORMATION

Follow all safety and processing instructions of the adhesive manufacturer.

### 10V8:

- 1) Fabricate a laminated forearm socket using the matching lamination form (here  $\varnothing$  50 mm).
- 2) Loosen the 4 attachment screws of the lamination ring.
- 3) Determine the position of the ratchet type rotation wrist unit in the socket so that it will be possible to comfortably actuate the levers with the contralateral side. Transfer the positions of the 4 bore holes of the lamination ring to the forearm socket (Fig. 3).
- 4) Glue the lamination ring into the forearm socket using 636K18 Sealing Resin Compact Adhesive and let cure. Then mount the ratchet type rotation wrist unit.

For mounting the 10V26 Flexion Device, the 4 screws on the 10V8 ratchet type rotation wrist unit used for attachment of the lamination ring must be disassembled. The lamination ring must be removed. Insert the wrist unit into the flexion device and fasten with the four screws (Fig. 4).

### 10V30:

- 1) Fabricate a laminated forearm socket using the matching lamination form (here  $\varnothing$  50 mm).
- 2) Determine the position of the ratchet type rotation wrist unit in the socket so that it will be possible to comfortably actuate the levers with the contralateral side.
- 3) Recess a slot for the release lever into the socket (Fig. 5). Unscrew the release lever for this purpose.

- 4) Glue the ratchet type rotation wrist unit into the forearm socket and let cure. Make sure that no adhesive enters the slot and the mechanism behind the slot.

## **5 Disposal**

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

## **6 Legal information**

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### ***6.1 Liability***

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### ***6.2 CE conformity***

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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