8E70=\*, 8E71=\*, 8E72=\*

<b>DE</b> Gebrauchsanweisung (Benutzer) .....	3
<b>EN</b> Instructions for use (user) .....	29
<b>FR</b> Instructions d'utilisation (Utilisateur) .....	55
<b>IT</b> Istruzioni per l'uso (Utilizzatore) .....	82
<b>ES</b> Instrucciones de uso (Usuario) .....	109
<b>PT</b> Manual de utilização (Usuário) .....	135
<b>NL</b> Gebruiksaanwijzing (Gebruiker) .....	162
<b>SV</b> Bruksanvisning (Användare) .....	189
<b>DA</b> Brugsanvisning (Bruger) .....	214
<b>NO</b> Bruksanvisning (Bruker) .....	239
<b>FI</b> Käyttöohje (Käyttäjä) .....	264
<b>HU</b> Használati utasítás (felhasználó) .....	290
<b>TR</b> Kullanma talimatı (Kullanıcı) .....	315
<b>RU</b> Руководство по применению (Пользователь) .....	340

**INFORMATION**

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## INFORMATION

Date of last update: 2022-06-30

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.


The products "8E70=\* bebionic hand EQD, 8E71=\* bebionic hand Short Wrist and 8E72=\* bebionic hand Flex" are referred to as the product/terminal device/hand in the following. These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

## 2 Product description

### 2.1 Design

The product consists of the following components:

1. Distal phalanx
2. Middle finger joint
3. Program switch (see page 33)
4. Back of the hand
5. Wrist joint
6. Distal thumb joint

### 2.2 Function

The product is a myoelectrically controlled and multiarticulating hand prosthesis.

The product enables objects to be gripped, pushed, pulled and carried with multiarticulating functionality. The product is designed to resemble the anatomical structure and weight of a human hand.

The thumb can be set to two different positions (opposition and lateral), allowing for 14 different grips (see page 29). The O&P professional can pre-configure 8 grips using the adjustment software.

A finger chain that protects the drive serves as overload protection for the 4 actively powered fingers. In case of overloading, the corresponding finger cannot be flexed anymore since the connection to the drive was severed. The finger chain can be replaced by an authorised O&P professional.

### 2.2.1 Grips

#### Moving the thumb from the lateral to the opposition position

- 1) Firmly grasp the base of the thumb with the free hand.
- 2) Apply controlled pressure to push the thumb inwards until it is opposite the palm.

#### Moving the thumb from the opposition to the lateral position

- 1) Firmly grasp the base of the thumb with the free hand.
- 2) Apply controlled pressure to push the thumb outwards until it is lateral to the palm.

#### Grips with opposed thumb (the thumb is opposite the palm)

##### Tripod grip

The index finger and middle finger are closed simultaneously with the thumb until the three fingers touch. Ring and little fingers continue to close until they meet resistance or the close signal stops. The thumb position has to be manually adjusted by the O&P professional for this grip so it can be carried out precisely.

Application examples: this grip can be used to grasp and hold objects (e.g. pens, coins).

##### Power grip

All fingers close until they encounter an object or the close signal stops. The thumb then moves towards the palm as well.

Application examples: this grip can be used to hold round objects (e.g. fruit, balls, glasses).

### **Finger adduction**

By spreading the fingers, a flat, thin object (< 3 mm/< 0.12 inch) can be held between the finger joints when closing the hand. Spreading the fingers is most effective with the power grip. It can also be used with the key grip and finger point. This grip **cannot** be selected in the adjustment software. Application examples: makes it possible to hold thin objects effectively (e.g. magazine, cutlery, toothbrush).

### **Hook grip**

This grip corresponds to a power grip with the fingers only partly closed. It makes it possible to engage objects with carrying loops and can also be initiated from the relaxed hand position. This grip **cannot** be selected in the adjustment software. Application examples: this grip makes it possible to carry bags.

### **Active index grip**

All fingers close and can pick up an object while the index finger remains extended. Subsequently the index finger can be individually flexed or extended by the user. An open signal causes the index finger to extend first followed by the remaining fingers. The object is released from the hand.

Application examples: this grip can be used for example to operate spray bottles.

### **Pinch grip**

For this grip, only the index finger and thumb meet while the remaining fingers close.

The thumb has to be manually adjusted by the O&P professional for this grip so it can be carried out precisely.

Application examples: this grip pattern makes it possible to grasp small objects (e.g. house keys, coins, closures, pens).

### **Precision closed grip**

The middle finger, ring finger and little finger are closed. The thumb moves to a half-closed position. Subsequently the index finger can be individually flexed or extended.

The thumb position has to be manually adjusted by the O&P professional for this grip so it can be carried out precisely.

Application examples: this grip makes it possible to grasp small objects and makes working at a table easier.

### **Precision open grip**

The middle finger, ring finger and little finger remain open. The thumb moves to a half-closed position. Subsequently the index finger can be individually flexed or extended. The thumb has to be manually adjusted by the O&P professional for this grip so it can be carried out precisely.

Application examples: this grip makes it possible to grasp small objects.

## Grips with the thumb in the lateral position (the thumb is sideways to the palm)

### Column grip

The thumb closes towards the palm, then the remaining fingers bend over the thumb. In doing so, the thumb blocks the closing movement of the index finger and lets it protrude beyond the other three fingers.

Application examples: this grip can be used to push large buttons (e.g. light switches) and for putting on clothes.

### Finger point

The middle finger, ring finger, little finger and thumb move towards the palm. The index finger remains extended.

Application examples: small buttons can be operated with this grip (e.g. keyboard, remote control).

### Key grip

The fingers close partway. In doing so, the thumb grasps the index finger sideways. This makes it possible to secure and release flat objects with the thumb, without movement of the remaining fingers.

Application examples: this grip can be used to hold thin objects without the remaining fingers moving (e.g. a spoon, paper, plate, credit card or key).

### Mouse grip

The thumb and little finger close to grasp the computer mouse sideways. Only when the thumb encounters resistance can the index finger be flexed. The index finger extends automatically when no close signal is applied. An open signal loosens the grip.

Application examples: this grip makes it possible to operate a computer mouse.

### Open palm grip

In the open hand position, the thumb is in the lateral position to achieve the largest possible flat palm. This grip **cannot** be selected in the adjustment software.

Application examples: this grip makes it possible to carry plates on the flat hand.

### Relaxed hand position

The thumb is in the lateral position and closes partway towards the palm. All fingers assume a slightly flexed position. When the close signal is repeated, the fingers move to form the hook grip. Application examples: this grip is recommended when the hand is not being actively used.

### 2.2.2 Program switch

A program switch is located on the back of the product (item 4). Both the program switch and flexion are not visible when using a prosthetic glove and have to be felt.

The switch incorporates various functions:

- Switching the terminal device on and off (see page 44)
- Switching between primary and secondary grips (see page 44)
- Activating donning mode (see page 45)

Depending on whether the terminal device is switched on or off and how long the program switch is pressed, the following functions can be carried out:

#### Terminal device switched on

Duration of pressing	Function	Beep signal	Vibration signal
Approx. 1 second	Switching between primary and secondary grips	1x short after releasing the program switch	1x short after releasing the program switch
Between 2 and 3 seconds	Switching the hand off	–	–
Longer than 4 seconds	Switching off the Bluetooth function	1x short	1x short
Longer than 4 seconds	Switching on the Bluetooth function	2x long	2x long

#### Terminal device switched off

Duration of pressing	Function	Beep signal	Vibration signal
Between 2 and 3 seconds	Switching the hand on	1x short after releasing the program switch	1x short after releasing the program switch
Approx. 3 seconds (until thumb opens)	Deactivating donning mode	–	–
Approx. 5 seconds (until thumb closes)	Activating donning mode	1x short	1x short

### 2.2.3 Switching modes

Switching between the default and alternative grip can be carried out as follows, depending on the selected mode:

- Mode 0: program switch
- Modes 1 through 4: with another open signal after fully opening the hand
- Mode 5: co-contraction signal after fully opening the hand

### 2.2.4 Wrist joint versions

The products "8E70=\*, 8E71=\*, 8E72=\* bebionic hand" are differentiated by different versions of the wrist joint:

### **8E70=\* bebionic hand EQD (with quick-disconnect wrist unit)**

Permits easy separation of the terminal device from the prosthetic socket. The terminal device can be removed quickly when needed via a 360° rotating movement and replaced by a different terminal device with the same quick-disconnect unit.

### **8E71=\* bebionic hand Short Wrist**

Low-profile connection for users with a long forearm or transcarpal amputation. The hand can be rotated against a constant friction, which can be adjusted during the fitting process.

### **8E72=\* bebionic hand Flex**

The flexion joint with quick-disconnect wrist unit makes positioning possible for the user at 20° or 40° of flexion, in the neutral position and at 20° or 40° of extension. The terminal device can be removed quickly when needed via a 360° rotating movement and replaced by a different terminal device with the same quick-disconnect unit.

## **3 Intended use**

### **3.1 Indications for use**

The product is intended **exclusively** for upper limb exoprosthetic fittings.

### **3.2 Conditions of use**

The product is intended **exclusively** for use on **one** user. Use of the product by another person is not approved by the manufacturer.

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports which involve excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing,

paragliding, etc.). Furthermore, the product should not be used to operate motor vehicles, heavy equipment (e.g. construction machines), industrial machines or motor-driven equipment. Permissible ambient conditions are described in the technical data (see page 48).

### 3.3 Indications

- Amputation level below-elbow, above-elbow and shoulder disarticulation
- For unilateral or bilateral amputation
- Dismelia of the forearm or upper arm
- The user must be able to understand usage instructions and safety notices and put them into practice.
- The user must fulfil the physical and mental requirements for perceiving optical/acoustic signals and/or mechanical vibrations

### 3.4 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

### 3.5 Qualification

The fitting of a patient with the product may only be carried out by O&P professionals who have been authorised with the corresponding Ottobock training.

## 4 Safety

### 4.1 Explanation of warning symbols

 <b>WARNING</b>	Warning regarding possible serious risks of accident or injury.
 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

### 4.2 Structure of the safety instructions

 <b>WARNING</b>
<b>The heading describes the source and/or the type of hazard</b>
The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:
> E.g.: Consequence 1 in the event of failure to observe the hazard
> E.g.: Consequence 2 in the event of failure to observe the hazard
▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

### 4.3 General safety instructions

 <b>WARNING</b>
<b>Non-observance of safety notices</b>
Personal injury/damage to the product due to using the product in certain situations.
▶ Observe the safety notices and the stated precautions in this accompanying document.

 <b>WARNING</b>
<b>Operating a vehicle with the prosthesis</b>
Accidents due to unexpected actions of the prosthesis.
▶ The prosthesis should not be used for the operation of motor vehicles and heavy equipment (e.g. construction machines).

**⚠ WARNING**

**Using the prosthesis while operating machines**

Injury due to unexpected actions of the prosthesis.

- ▶ The prosthesis should not be used to operate industrial machines or motor-driven equipment.

**⚠ WARNING**

**Use of damaged power supply unit, adapter plug or battery charger**

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

**⚠ WARNING**

**Use of the product while handling a firearm.**

Injury due to unexpected product behaviour.

- ▶ The product may not be used for handling firearms.

**⚠ WARNING**

**Skin contact with leaking lubricants due to mechanical defects**

Injury due to skin irritation.

- ▶ Avoid contact between leaking lubricants and the mouth, nose and eyes.
- ▶ The product must be inspected by an authorised Ottobock Service Center.

**⚠ WARNING**

**Operating the prosthesis near active implanted systems**

Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic emissions of the prosthesis.

- ▶ When operating the prosthesis in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

**⚠ CAUTION**

**Signs of wear on the product**

Injury due to faulty control or malfunction of the product

- ▶ In the interest of personal safety and in order to maintain operating reliability, the prosthesis has to be inspected by a O&P professional in case of noticeable restrictions of functionality.
- ▶ Note that the functionality of the terminal device may be limited if the battery charge level gets too low.

**⚠ CAUTION**

**Use of a damaged product**

Injury due to loss of product functionality.

- ▶ Prior to use, conduct an external visual inspection to verify that all parts of the product are undamaged.

- ▶ In case of damage, have the product repaired promptly.

### **⚠ CAUTION**

#### **Prosthesis use when battery charge level is too low**

Injury due to unexpected prosthesis behaviour

- ▶ Check the current charge level before use and charge the prosthesis if required.
- ▶ Note that the operating time of the prosthesis may be reduced at low ambient temperatures or due to ageing of the battery.
- ▶ Note that the actions/reactions of the terminal device become slower when the battery voltage gets very low.
- ▶ Note that only a few grips or actions will still be possible with the terminal device when the battery voltage gets very low.
- ▶ A small opening width may indicate a low battery voltage.

### **⚠ CAUTION**

#### **Penetration of product with dirt and humidity**

Injury due to unexpected product behaviour or malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

### **⚠ CAUTION**

#### **Independent manipulation of the product or product components**

Injury due to faulty control or malfunction of the product due to manipulation

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The battery may only be handled by authorised, qualified Ottobock personnel (no replacement by the user).
- ▶ Opening and repairing the product and any damaged components, and removing the back of the hand, is reserved exclusively for authorised, qualified Ottobock personnel.

### **⚠ CAUTION**

#### **Use of unapproved components**

> Injury due to product malfunction as a result of reduced interference resistance.

> Interference of other electronic devices due to increased emissions.

- ▶ Use the product only in combination with the components, signal converters and cables listed in the section "Scope of delivery" (see page 41).

### **⚠ CAUTION**

#### **Failure to observe the instructions for use of all prosthesis components used**

Injury due to unexpected product behaviour.

- ▶ Observe the instructions for use of the prosthesis components used.

### **⚠ CAUTION**

#### **Insufficient skin contact of the electrodes**

Injury due to unexpected product behaviour as a result of uncontrolled activation of the prosthesis component.

- ▶ The electrodes are to be placed on intact skin only and with as much electrode-skin contact as possible.
- ▶ Make sure that skin contact of the electrodes continues to be given when carrying heavy loads.

- ▶ If the product cannot be controlled correctly using muscle signals, turn off the entire prosthesis and contact your O&P professional.

#### 4.4 Information on the Power Supply/Battery Charging

##### **⚠ WARNING**

##### **Charging the prosthesis without taking it off**

Risk of electric shock due to defects in the power supply unit or in the battery charger.

- ▶ For safety reasons, remove the prosthesis prior to charging it.

##### **⚠ CAUTION**

##### **Charging the product with damaged power supply unit/charger/charger cable**

Injury due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Check the power supply unit, charger and charger cable for damage before use.
- ▶ Replace any damaged power supply unit, charger or charger cable.

##### **NOTICE**

##### **Use of incorrect power supply unit/battery charger**

Damage to product due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for this product by Ottobock (see instructions for use and catalogues).

#### 4.5 Battery charger information

##### **NOTICE**

##### **Penetration of dirt and humidity into the product**

Lack of proper charging functionality due to malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

##### **NOTICE**

##### **Mechanical stress on the power supply/battery charger**

Lack of proper charging functionality due to malfunction.

- ▶ Do not subject the power supply/battery charger to mechanical vibrations or impacts.
- ▶ Check the power supply/battery charger for visible damage before each use.

##### **NOTICE**

##### **Operating the power supply unit/charger outside of the permissible temperature range**

Lack of proper charging functionality due to malfunction.

- ▶ Only use the power supply unit/charger for charging within the allowable temperature range. The section "Technical data" contains information on the allowable temperature range (see page 48).

##### **NOTICE**

##### **Independent changes or modifications carried out to the battery charger**

Lack of proper charging functionality due to malfunction.

- ▶ Have any changes or modifications carried out only by Ottobock authorised, qualified personnel.

## 4.6 Information on Proximity to Certain Areas

### CAUTION

#### **Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm from HF communication devices is recommended.

### CAUTION

#### **Operating the product in very close proximity to other electronic devices**

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Do not operate the product in the immediate vicinity of other electronic devices.
- ▶ Do not stack the product with other electronic devices during operation.
- ▶ If simultaneous operation cannot be avoided, monitor the product and verify proper use in the existing setup.

### CAUTION

#### **Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

## 4.7 Information on Use

### CAUTION

#### **Mechanical stress on the product**

Injury due to faulty control or malfunction of the product.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

### CAUTION

#### **Improper handling**

Injury due to faulty operation or malfunction of the product.

- ▶ Seek expert guidance in the proper use of the product.

### CAUTION

#### **Improper product care**

- > Injuries due to faulty control/malfunction of the product or damage to the mechanical components

- > Damage or breakage due to embrittlement of plastics caused by the use of acetone, petrol or similar solvents.
- ▶ Clean the product only as described in the section "Cleaning and care" (see page 46).
- ▶ Do not clean the product under running water.
- ▶ When using a prosthetic glove, also note the instructions for use of the glove.

**⚠ CAUTION**

**Grasping objects with incorrect gripping forces**

Injury due to unexpected product behaviour.

- ▶ Note that the gripping has to be controlled manually depending on the consistency (soft/hard) of the object being grasped.

**⚠ CAUTION**

**Overloading due to unusual activities**

Injury due to unexpected product behaviour as a result of a malfunction.

- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures your personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately by an O&P professional. If necessary, he or she will forward the product to an authorised Ottobock Service Centre.

**⚠ CAUTION**

**Risk of pinching between the fingertips**

Injury due to pinching of body parts.

- ▶ Ensure that no body parts are between the fingertips when using the product.
- ▶ Ensure that no body parts are between the fingertips when closing the hand.
- ▶ When closing the hand, ensure that fingers and other body parts are not in the area of the finger joints.
- ▶ Make sure the product is switched off for cleaning.

**⚠ CAUTION**

**Insufficient distance from sources of high heat**

Inflammation of the product.

- ▶ Do not expose the product to sources of high heat (fire, stovetop, space heater, radiator etc.).
- ▶ Do not grasp or hold any red-hot objects with the product.

**⚠ CAUTION**

**Unintentional unlocking of the terminal device**

Injury due to releasing the terminal device from the forearm (e.g. while carrying objects).

- ▶ When connecting the hand to the prosthetic socket or components, ensure that the connection is carried out correctly.

## INFORMATION

Avoid exposing the product without a prosthetic glove to direct sunlight or UV light (solarium) for extended periods.

## 5 Scope of Delivery and Accessories

### 5.1 Scope of delivery

- 1 pc. cosmetic case for battery charger and power supply
- 1 pc. prosthesis passport
- 1 pc. Instructions for use (user)

#### If no 12K100N\* elbow component is used

- 1 pc. 757L35 MyoCharge Integral battery charger (including 757L16-4 power supply)

#### If the 12K100N\* elbow component is used

- 1 pc. 757L24 DynamicArm battery charger

#### Terminal devices

- 1 pc. 8E70=\* bebionic hand EQD
- or
- 1 pc. 8E71=\* bebionic hand Short Wrist
- or
- 1 pc. 8E72=\* bebionic hand Flex

### 5.2 Accessories

- 8S710=\* prosthetic glove
- 8S711=\* prosthetic glove

## 6 Use

### 6.1 Charging the prosthesis battery

## INFORMATION

#### Use of the 12K100N=\* elbow component

When using this elbow component, note the information on charging the battery, querying the charge level and the feedback (beeps and vibration signals) in the instructions for use (user) included with this elbow component.

#### 6.1.1 Charging the battery

## ⚠ WARNING

#### Charging the prosthesis without taking it off

Risk of electric shock due to defects in the power supply unit or in the battery charger.

- ▶ For safety reasons, remove the prosthesis prior to charging it.

The following points must be observed when charging the battery:

- The capacity of a fully charged battery is sufficient for one full day.
- We recommend charging the product every day when used on a daily basis.
- The battery should be charged for at least 3 hours prior to initial use.

### 6.1.1.1 Connect battery charger to the product

- 1) Connect the charging plug to the charging receptacle on the prosthetic socket.  
→ A correct connection between the battery charger and the prosthesis is indicated by feedback (Status signals).
- 2) The charging process starts.  
→ The prosthesis is switched off automatically.
- 3) Disconnect the prosthesis when the charging process is complete.

### 6.1.2 Safety Shutoff

The purpose of the prosthesis safety shutoff is to protect the battery; it is triggered in case of:

- Excessively high or low temperature
- Overvoltage or undervoltage
- Short circuit

After a short circuit, the charging plug must be plugged into the charging receptacle and then unplugged in order for the electronics to be activated.

### 6.1.3 Display of the current charge level

The battery charge level can be queried at any time.

- 1) With the prosthesis switched on, press the charging receptacle button and hold for less than one second.
- 2) The LED display on the charging receptacle provides information on the current battery charge level.

Charging receptacle	Event
●	Battery fully charged (illuminated in green)
●	Battery 50% charged (illuminated in yellow)
●	Battery drained (illuminated in orange)

### 6.1.4 Beep signals

Beep signal	Additional display	Event
1 x long	-	<ul style="list-style-type: none"> <li>• Turn off prosthesis on charging receptacle</li> <li>• Charging begins (charging plug connected to charging receptacle)</li> <li>• Charging ends (charging plug disconnected from charging receptacle)</li> </ul>
2 x short	LED on charging receptacle lights up briefly	Switch on prosthesis
3 x short	-	Battery voltage too low, prosthesis shuts off automatically

### 6.1.5 Behaviour when the battery is low

If the product battery is too low, the speed of prosthesis control becomes slower. Executing a certain grip takes longer than usual, and ultimately the hand switches off entirely and remains in the open state.

## INFORMATION

When the grips slow down, objects held in the hand have to be released since the terminal device ultimately switches off and then cannot execute grips anymore.

### 6.2 Thumb adjustment

#### Moving the thumb from the lateral to the opposition position

- 1) Firmly grasp the base of the thumb with the free hand.
- 2) Apply controlled pressure to push the thumb inwards until it is opposite the palm.

#### Moving the thumb from the opposition to the lateral position

- 1) Firmly grasp the base of the thumb with the free hand.
- 2) Apply controlled pressure to push the thumb outwards until it is lateral to the palm.

### 6.3 Unlocking/locking the wrist joint (8E72=\*)

The individual flexion and extension of the wrist joint can be locked in five different positions (in 20° increments each).

- 1) Press the release button in the direction of the arrow.
- 2) Move the terminal device to the desired position while holding the release button. From the relaxed hand position, it engages at 20° and 40° in each direction.
- 3) Let go of the release button to lock the terminal device in the respective position.

## 6.4 Applying/removing the terminal device

### INFORMATION

Applying/removing the terminal device works only with the terminal devices "8E70=\* bebionic hand EQD" and "8E72=\* bebionic hand Flex".

## Separating the terminal device from the socket

### ⚠ CAUTION

#### Unintentional unlocking of the terminal device

Injury due to releasing the terminal device from the forearm (e.g. while carrying objects).

► When connecting the hand to the prosthetic socket or components, ensure that the connection is carried out correctly.

- 1) Turn the prosthesis off with the button in the charging receptacle (prosthetic socket) or the switch on the elbow joint.
- 2) Rotate the terminal device once around its own axis until a slight resistance is felt (about 360°).
- 3) Overcome this resistance and pull the terminal device off the socket.

#### Attaching the terminal device to the socket

- 1) Insert the quick-disconnect wrist unit into the lamination ring and press it in firmly.
- 2) Turn the terminal device slightly to the left or right.
- 3) Verify proper attachment by pulling on the terminal device.

### INFORMATION

Please see the 646T332 technical information for proper handling with the lamination ring.

## 6.5 Terminal device on/off

### Switching the hand on

- 1) Press and hold the program switch on the back of the hand for no longer than 2 to 3 seconds.
- 2) After releasing the program switch, a single short feedback signal is generated (see page 51).

**INFORMATION: If a feedback signal was already generated in step 1, the Bluetooth function has been activated or deactivated.**

→ The terminal device is switched on.

After the prosthesis has been switched on with the button in the charging receptacle or the switch on the elbow joint, the terminal device is switched on as well.

### Switching the hand off

- Press and hold the program switch on the back of the hand for no longer than 2 to 3 seconds.
- If no feedback signal is generated after releasing the program switch, only the terminal device is switched off. Other prosthesis components, such as an elbow joint or electric rotator, can still be used.

**INFORMATION: If a feedback signal is generated after releasing or while pressing the program switch, or if the thumb closes towards the palm, the program switch was not pressed long enough or pressed too long.**

When the prosthesis is switched off with the button in the charging receptacle or the switch on the elbow joint, all other prosthesis components including the terminal device are switched off.

## 6.6 Switching between primary and secondary grips

- Briefly press the program switch. A short single feedback signal is generated after it is released (see page 51).
- The product has switched between the primary and secondary grips or vice versa.

## 6.7 Switching between default and alternative grips

Switching between the default and alternative grip can be carried out as follows, depending on the selected mode:

- Mode 0: program switch
- Modes 1 through 4: with another open signal after fully opening the hand
- Mode 5: co-contraction signal after fully opening the hand

## 6.8 Donning mode

### Activating donning mode

- 1) Move the thumb to the opposed position.
- 2) With the terminal device switched off, press and hold the program switch on the back of the hand until the thumb automatically moves to donning mode.  
→ The thumb moves in towards the palm.  
or
- 1) Switch the prosthesis on with the button in the charging receptacle or the switch on the elbow joint.
- 2) During the initialisation phase of the terminal device, press and hold the program switch on the back of the hand until the thumb automatically moves to donning mode.  
→ The thumb moves in towards the palm.

### INFORMATION

#### **Sending the product to an authorised Ottobock Service Centre**


Set the thumb to the lateral position and then activate "donning mode" on the terminal device.

### Deactivating donning mode

- ▶ With the thumb pointing inward, press and hold the program switch on the back of the hand until the thumb opens.
- Donning mode is deactivated and the terminal device is switched on.

## 7 User information

### 7.1 Hand Rotation

## 7.2 Supporting the hand

### INFORMATION

Supporting oneself with the product "8E72=\* bebionic hand Flex" may lead to buckling of the wrist joint.

## 7.3 Using Finger Adduction/Abduction, e. g. when Carrying a Bag

## 8 Cleaning and care

- 1) Switch off the product before cleaning.
- 2) Clean the product with a damp cloth and mild soap when needed.  
Make sure that no liquids get into the product and product components.
- 3) Dry the product with a lint-free cloth and allow it to air dry fully.

### INFORMATION

The hand can be cleaned three times per day on average.

### INFORMATION

When using a prosthetic glove, note the cleaning information in the instructions for use for the prosthetic glove.

## 9 Legal information

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for

damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

## 9.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

## 9.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

This product meets the requirements of the 2014/53/EU directive.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

## 9.4 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use in this chapter.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and
- 2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

### **Caution: Exposure to Radio Frequency Radiation.**

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s)..

Operation is subject to the following two conditions:

- (1) This device may not cause interference.

(2) This device must accept any interference, including interference that may cause undesired operation of the device.

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**Caution: Exposure to Radio Frequency Radiation.**

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population.

**Caution: Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.**

**10 Technical data**

<b>Ambient conditions</b>	
Storage in original packaging	+5°C/+41°F to +40°C/+104°F Max. 85% relative humidity, non-condensing
Transport in original packaging	-25°C/-13°F to +70°C/+158°F Max. 90% relative humidity, non-condensing
Storage and transport without packaging	-25°C/-13°F to +70°C/+158°F Max. 90% relative humidity, non-condensing
Operation	-5°C/+23°F to +45°C/+113°F Max. 95% relative humidity, non-condensing

<b>General information</b>	<b>bebionic hand EQD</b>	<b>bebionic hand Short Wrist</b>	<b>bebionic hand Flex Wrist</b>
Reference number	8E70=*	8E71=*	8E72=*
Weight of the bebionic hand small	433 g/0.95 lbs	Approx. 402 g/0.89 lbs	Approx. 504 g/1.1 lbs
Weight of the bebionic hand medium	616 g/1.36 lbs	Approx. 588 g/1.3 lbs	Approx. 689 g/1.52 lbs
Opening width (between index finger and opposed thumb)	75 mm		
Extension/flexion of the wrist	-	-	-40° to +40° in 20° increments
Service life	5 years		
Version of the terminal device	Hardware and firmware version accessible to the O&P professional in the adjustment software		

<b>Load limits</b>	
Force on single finger (static)	32 N
Transverse force on single finger (static)	44 N
Force on chassis (static, supporting the hand)	500 N
Force with closed hand (static, carrying a bag)	152 N

<b>Load limits</b>	
Forces on thumb (static)	40 N

<b>Data transfer</b>	
Wireless technology	Bluetooth Smart/Low Energy
Range	min. 2 m / 6.7 ft
Frequency range	2402 MHz to 2480 MHz
Modulation	GFSK
Maximum output power	9.6 dBm

**If no 12K100N\* elbow component is used**

<b>Rechargeable battery of the product</b>	<b>757B35=3</b>	<b>757B35=4</b>	<b>757B35=5</b>
Temperature range for charging [°C]	+5 °C to +40 °C		
Temperature range for charging [°F]	+41 °F to 104 °F		
Capacity [mAh]	1150	2350	3450
Rechargeable battery lifetime [years]	2		
Product behaviour during the charging process	The product is non-functional		
Operating time of the product with fully charged battery during average use [hours]	Approx. 8	Approx. 16	Approx. 24
Charging times (rechargeable battery fully discharged) [hours]	Approx. 2.5	Approx. 3	Approx. 3
Nominal voltage [V]	Approx. 7.4		
Battery technology	Li-ion		
Approved battery chargers	757L35		

<b>Battery charger</b>	
Reference number	757L35
Product service life	5 years
Max. charging current	1 A
Max. charging voltage	10 V
Battery charger operating voltage	12 V
Power supply operating voltage	100 V - 240 V
Power supply operating frequency	50 Hz - 60 Hz

<b>Power supply unit</b>	
Reference number	757L16-4
Type	FW8001M/12
Storage and transport in original packaging	-40 °C/-40 °F to +70 °C/+158 °F 10% to 95% relative humidity, non-condensing
Storage and transport without packaging	-40 °C/-40 °F to +70 °C/+158 °F 10% to 95% relative humidity, non-condensing




<b>Power supply unit</b>	
Operation	0 °C/+32 °F to +50 °C/+122 °F Max. 95% relative humidity Air pressure: 70–106 kPa (up to 3,000 m without pressure equalisation)
Input voltage	100 V~ to 240 V~
Mains frequency	50 Hz to 60 Hz
Output voltage	12 V $\overline{=}$

**If the 12K100N\* elbow component is used**

<b>Battery charger</b>	
Reference number	757L24
Type	2240
Product service life	5 years
Max. charging current	1.2 A
Max. charging voltage	12.6 V
Operating voltage	100 V - 240 V~
Operating frequency	50 Hz - 60 Hz

## 11 Appendix

### 11.1 Product dimensions

		<b>bebionic hand small</b>	<b>bebionic hand medium</b>
Palm of the hand		72 mm	85 mm
Max. hand width		122 mm	136 mm
Hand length including fingers		162 mm	188 mm
Finger length		75 mm	91 mm

## 11.2 Symbols Used

Manufacturer

Type BF applied part

Compliance with the requirements according to "FCC Part 15" (USA)

Compliance with the requirements under the "Radiocommunications Act" (AUS)

Non-ionising radiation

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

Serial number (YYYY WW NNN)  
YYYY – year of manufacture  
WW – week of manufacture  
NNN – sequential number

Article number

Medical device

Protect from moisture

## 11.3 Operating states/error signals

The prosthesis indicates operating states and error messages through beeps and vibration signals.

## INFORMATION

### Use of the 12K100N=\* elbow component

When using this elbow component, note the information on charging the battery, querying the charge level and the feedback (beeps and vibration signals) in the instructions for use (user) included with this elbow component.

#### 11.3.1 Beep and vibration signals

## INFORMATION

### Feedback signals can be turned off

When the feedback signals have been turned off in the adjustment software, the beeps and/or vibration signals are not generated in certain cases (see table). Signals in case of a product error are generated even when the feedback signals have been turned off.

Beep signal	Vibration signal	When	Signal can be deactivated	Function
1x short	1x short	After releasing the program switch	Yes	Switching between primary and secondary grips
1x short	1x short	After releasing the program switch	Yes	Hand is switched on
1x short	1x short	While pressing the program switch	No	Bluetooth function is deactivated
1x short	1x short	While pressing the program switch	Yes	Donning mode was activated
2x long	2x long	While pressing the program switch	No	Bluetooth function is activated
Long for 3 seconds	Long for 3 seconds	During product initialisation	No	Error, the product must be inspected by a prosthetist.

#### 11.4 Directives and manufacturer's declaration

##### 11.4.1 Electromagnetic environment

This product is designed for operation in the following electromagnetic environments:

- Operation in a professional healthcare facility (e.g. hospital, etc.)
- Operation in areas of home healthcare (e.g. use at home, use outdoors)

Observe the safety notices in the section "Information on proximity to certain areas" (see page 39).

**Table 1 – Electromagnetic emissions for all devices and systems**

Interference measurements	Compliance	Electromagnetic environment directive
---------------------------	------------	---------------------------------------

HF emissions according to CISPR 11	Group 1	The product uses HF energy exclusively for its internal functioning. Its HF emissions are therefore very low, and interference with neighbouring electronic devices is unlikely.
HF emissions according to CISPR 11	Class B	The product is intended for use in other than residential facilities and in facilities connected directly to a public low-voltage network that also supplies residential buildings.
Harmonics according to IEC 61000-3-2	N/A	
Voltage fluctuations/flicker according to IEC 61000-3-3	N/A	

**Table 4 – casing**

Phenomenon	EMC basic standard or test procedure	Interference immunity test level	
		Professional healthcare facility	Environments in areas of home healthcare *)
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air,	
High frequency electromagnetic fields	IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz	12 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz
High frequency electromagnetic fields in the immediate vicinity of wireless communication devices	IEC 61000-4-3	See Table 9	
Magnetic fields with rated power frequencies	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	

\*) Tests conducted

**Table 5 – alternating current gate for the supply input**

Phenomenon	EMC basic standard or test procedure	Interference immunity test level	
		Professional healthcare facility	Environments in areas of home healthcare
Electrical fast transients/bursts	IEC 61000-4-4	± 2 kV 100 kHz repetition rate	
Surges Line against line	IEC 61000-4-5	± 0.5 kV, ± 1 kV	
Conducted interference induced by high-frequency fields	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM frequency bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz	3 V/ 0.15 MHz to 80 MHz 6 V in ISM and amateur radio frequency bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz

Voltage drops	IEC 61000-4-11	0 % $U_T$ ; 1/2 period at 0, 45, 90, 135, 180, 225, 270 and 315 degrees
		0 % $U_T$ ; 1 period and 70 % $U_T$ ; 25/30 periods Single phase: at 0 degrees
Voltage interruptions	IEC 61000-4-11	0 % $U_T$ ; 250/300 periods

**Table 9 – Test specifications for the interference resistance of casings against high frequency wireless communication devices**

Test frequency [MHz]	Frequency band [MHz]	Radio service	Modulation	Maximum power [W]	Distance [m]	Interference immunity test level [V/m]
385	380 to 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 to 470	GMRS 460, FRS 460	FM $\pm 5$ kHz deviation 1 kHz sine	1.8	0.3	28
710	704 to 787	LTE band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 to 960	GSM 800/90- 0, TETRA 800, iDEN 820, CDMA 850, GSM 800/90- 0, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1720	1700 to 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1845						
1970						
2450	2400 to 2570	Bluetooth WLAN 802.1- 1 b/g/n, RFID 2450 LTE band 7	Pulse modulation 217 Hz	2	0.3	28
5240	5100 to 5800	WLAN 802.1- 1 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5500						
5785						

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- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

— Consult the dealer or an experienced radio/ TV technician for help.

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This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s).

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(1) This device may not cause interference.

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