

28U22 WalkOn Flex Junior

Active dynamic gait and stability for children

Features & Benefits

- Improves stability for children with drop foot
- Provides a more physiological and symmetrical gait with fluid rollover
- Dynamic design and advanced carbon fiber construction for excellent energy storage and return during gait cycle
- Medial guidance of the longitudinal arch prevents supination in solid footwear
- Spiral design provides greater flexibility at heel strike
- Lightweight, low-profile, durable
- Trimmable footplates easily shaped with scissors, often requires only one office visit
- New and improved calf band with phase change material (temperature regulative)
- Can be used with other dynamic ankle-foot orthoses (DAFOs)
- No activity or weight restrictions
- 60-day money back guarantee and 12-month warranty (tested for one year durability)
- Packaged in a bright, kid-friendly sports bag


Indications

Light to moderate drop foot

Resulting from:

- Peroneal paralysis
- Stroke
- Traumatic brain injury
- Multiple sclerosis
- Neural muscular atrophy

Clinical indications must be determined by the practitioner.



Technical Specifications

Material Number	Children's Shoe Size	Side	Height
28U22=L24-27	7 - 9.5	Left	20.5 cm
28U22=R24-27	7 - 9.5	Right	20.5 cm
28U22=L27-30	9.5 - 12	Left	23.5 cm
28U22=R27-30	9.5 - 12	Right	23.5 cm
28U22=L30-33	12 - 14.5	Left	26.5 cm
28U22=R30-33	12 - 14.5	Right	26.5 cm
28U22=L33-36	3.5 - 6	Left	29.5 cm
28U22=R33-36	3.5 - 6	Right	29.5 cm